



**Winco**  
Medical Furnishings





## NON-TRENDELENBURG RECLINERS

5000 Series



**Winco**  
Medical Furnishings



# VERŌ

## TREATMENT RECLINER

PREMIUM SERIES



*Verō PRC with optional 3" dual casters, urethane arm caps and headrest cover & pillow set – shown in Bottle (85) and Limestone (57).*

Continuously evolving patient care requirements create opportunities for enhanced solutions. The Verō provides total flexibility in matching your application with exactly the features you require and the comfort your patients demand.

The Verō recliner's modular design allows clinicians to choose the best configuration to perfectly suit their specific needs:

- New Power Recline model (Verō PRC)!
- Multiple recline options and infinite position backrest promote patient independence.
- Bilateral and foot-operated controls provide clinician comfort and convenience.
- Durable construction and easy-to-clean surfaces enable increased efficiency.

**DESIGNED BY WINCO, BUILT BY YOU.  
THAT'S JUST SMART.**





# VERO CARE CLINER

PREMIUM SERIES

## ANATOMY OF THE VERO CARE CLINER

**Fastener-Free Headrest Cover**

*Optional for 5 Series*

**Push Handle**

*Option*

**Infinite Position Backrest**

**Power Recline**

*PRC Only*

**Gas Assist Control Lever**

*5X and 6X models only*

**Accessory Hook**

*Option*

**Power Recline Pendant**

*Option: PRC Only*

**Urethane Arm Cap**

*Option*

**Intuitive Side Tables**

*Option*

*Premium version shown*

**Foot Operated Swing Arm Release**

*Swing Arm Models only*

**Weight-Bearing, Retractable Foot Plate Option**

*5" Caster models only*

**Bilateral Trendelenburg Access**

*6 Series Models only*

**IV Pole + Holder**

*Option*

**Easily Accessed Heat & Massage or Heat Only Switch**

*Option*

50+ VINYL COLOR OPTIONS AVAILABLE

**Footrest with "Heel Drop" Comfort Zone**

**Independent Leg Rest**


**Durable, Easy-to-Clean, Powder-Coated Frame**

**180 Degree Swing Arm**

*Swing Arm models only*

**Fixed Front Casters with Accessible Lock Tabs**

*5" Casters only  
(Pedestal Feet or 3" Caster options available)*



## TRENDELENBURG

**Verō** Recliners are available with Trendelenburg (6 series) or without (5 series). On the **Verō**, Trendelenburg is controlled with bilateral access foot levers, reducing clinician's need to bend to reach controls. PRC models (with power recline) offer push-button recovery from Trendelenburg to minimize potential for clinician back strain.

Exclusive to **Verō**, Secure-Glide™ provides gentle patient positioning while avoiding injury to caregivers. The controlled descent helps eliminate any need or natural reflex to manually control the rapid patient descent in emergency situations.



Secure-Glide™ for stress-free positioning



Foot lever for exclusive Secure-Glide™ Trendelenburg descent


## RECLINE CONTROL OPTIONS

All **Verō** models feature convenient and comfortable infinite positioning of the back rest. With **Verō**, you have three options available for controlling recliner positioning:

**Power** Easy recline and recovery from trendelenburg at the push of a button. One switch standard. Dual switches and/or pendant options available.

**Gas Assist** provides positive locking position of the backrest angle with gas cylinder assisted operation. Locking action is controlled by easily accessed bilateral control levers.

**Push-Back** allows seamless positioning of the desired back rest angle. The gravity operated, gentle resistance designed into the chair maintains selected positioning.



**Power Recline**  
PRC models



**Gas Assist**  
5X and 6X models



**Push Back**  
5X and 6X models

## ARM OPTIONS

Functionality is enhanced through the availability of two arm configurations. Standard fixed arms are cost-effective, sturdy and durable. Optional swing arms allow 180 degree movement of the arm for convenient cleaning access and barrier-free patient transfer.

**Swing Arm Option** Swing arms are ideal when patient transfer and easy access are important considerations. Each side of the **Verō** is easily opened with convenient foot controls.

**Fixed Arm Option** Fixed arms are suitable for most general purpose operations when transfer of patient is not required or where working space is limited.



## ARM CAP OPTIONS

**Urethane Cap Option** Urethane arm caps provides a unique, upscale appearance and added durability.

**Urethane Arm Cap**  
Option





## MATCH YOUR PURPOSE

The **Verō** modular design allows choices in how the chair is anchored or moved. Whether stationary or for patient transport\*, **Verō** has a perfect solution.

**Pedestal Feet** When mobility is not required or desired, pedestal feet provide a less clinical look and a lower seat height.

**Casters** provide ease of movement for cleaning, maintenance and staging. Two size options are available. Both include lockable, foot lever release controls.

**3" Casters** Lower seated height with the flexibility of mobility.

**5" Casters** are the optimal choice for ease in movement.



*It is recommended that after moving the chair always lock the casters prior to patient ingress/egress.*

## FOOT PLATE

Patients will appreciate the added comfort of the retractable foot plate.

The optional footplate is available for models with 5" casters. Its sturdy design is weight-bearing.

The footplate locks as weight is applied to prevent unwanted or unexpected movement.

**Note: Patients should be assisted at all times by a staff member when using the footplate.**



*Footplate is easily positioned by foot.*

## PUSH HANDLE OPTION

The push handle option is recommended in combination with footplate models if the chair is to be used for in-facility patient transport\*.

*\*Only 5X, 6X or PRC models with footplate and push handle are suitable for patient transport. See User Manual for important details.*





## SIZE / CAPACITY OPTIONS FOR MAXIMIZED FLEXIBILITY

The **Verō** is available in two sizes: standard and X-large. No matter which size is specified, **Verō** offers comfort, beauty and durability.

- Standard seat width is 22.5" between arms with a 350 pound weight capacity. This chair makes petite to average size patients feel right at home.
- X-large seat width is 26.5" between arms with a 500 pound weight capacity, making it suitable to care for Bariatric patients and provide a comfortable experience.



X-Large 500 lb Capacity



Standard 350 lb Capacity

## TABLE OPTIONS

There are three table options; standard, pivot and premium. All tables on the **Verō** are available in left, right or both configurations and can be ordered in a variety of combinations to suit facility requirements and clinician preference. All feature a stain resistant, non-porous material and all but the pivot table have a lip edge design to contain small spills.



Standard Table



Pivot Table




Premium Table

Each standard table adds 2.75" (2.5" for XL) to the overall chair width when lowered. Each premium table adds 1.25" (1" for XL) to the overall chair width when lowered. This added width should be considered when ordering table options appropriate for your facility.

- **Standard table** This value priced lift-to-raise-and-lower table will be familiar to current Winco Elite series owners.
- **Pivot table** swivels so patient can have the table in front of them.
- **Premium table** features an intuitive single hand up/down feature with audible positive lock and release.

## COMPLEMENT DESIGN & UTILITY

- **Heat Only Option** provides three patient-selected soothing warmth ranges for maximum comfort.
- **Heat and Massage Option** adds a selection of five gentle vibrating massage modes to the heat option, enhancing patient comfort level.
- **Optional Accessory Hook** The stainless steel hook remains flush until needed, making it a convenient choice for a Foley bag. Available for either left, right or both sides.
- **Optional IV Pole and Holder** for either left, right or both sides.
- **Headrest Cover and Armrest Covers** Covers extend the life of the headrest and armrests and are fastener-free, making them easy to lift on and off.



**Heat & Massage**  
3 levels of heat, 5 levels of massage

**Heat Only**  
Low/Med/High

CONVENIENTLY  
PLACED CONTROLS  
FOR PATIENT  
INDEPENDENCE



**Accessory Hook**  
Shown here with Foley bag



**IV Holder Mount**  
Available for left, right or both sides

# MODELS & SPECIFICATIONS

MOBILITY OPTION	PEDESTAL FEET				3" CASTERS				5" CASTERS			
	STANDARD		X-LARGE		STANDARD		X-LARGE		STANDARD		X-LARGE	
	PUSH BACK/ GAS BACK	POWER	PUSH BACK/ GAS BACK	POWER	PUSH BACK/ GAS BACK	POWER	PUSH BACK/ GAS BACK	POWER	PUSH BACK/ GAS BACK	POWER	PUSH BACK/ GAS BACK	POWER
SEAT HEIGHT TO FLOOR	19.25 in (48.9 cm)	20.25 in (51.4 cm)	19.25 in (48.9 cm)	20.25 in (51.4 cm)	19.75 in (50.2 cm)	20.75 in (52.7 cm)	19.75 in (50.2 cm)	20.75 in (52.7 cm)	22.75 in (57.8 cm)	23.75 in (60.3 cm)	22.75 in (57.8 cm)	23.75 in (60.3 cm)
SEAT WIDTH	22.5 in (57.2 cm)	22.5 in (57.2 cm)	26.5 in (67.3 cm)	26.5 in (67.3 cm)	22.5 in (57.2 cm)	22.5 in (57.2 cm)	26.5 in (67.3 cm)	26.5 in (67.3 cm)	22.5 in (57.2 cm)	22.5 in (57.2 cm)	26.5 in (67.3 cm)	26.5 in (67.3 cm)
SEAT DEPTH	22 in (55.9 cm)	23.5 in (59.7 cm)	22 in (55.9 cm)	23.5 in (59.7 cm)	22 in (55.9 cm)	23.5 in (59.7 cm)	22 in (55.9 cm)	23.5 in (59.7 cm)	22 in (55.9 cm)	23.5 in (59.7 cm)	22 in (55.9 cm)	23.5 in (59.7 cm)
OVERALL HEIGHT	44.75 in (113.7 cm)	44 in (111.8 cm)	44.75 in (113.7 cm)	44 in (111.8 cm)	45.25 in (115 cm)	44.5 in (113 cm)	45.25 in (115 cm)	44.5 in (113 cm)	48.25 in (122.6 cm)	47.5 in (120.7 cm)	48.25 in (122.6 cm)	47.5 in (120.7 cm)
OVERALL WIDTH (FIXED ARMS / NO TABLES)*	29.5 in (11.61 cm)	29.5 in (11.61 cm)	34 in (13.39 cm)	34 in (13.39 cm)	29.5 in (11.61 cm)	29.5 in (11.61 cm)	34 in (13.39 cm)	34 in (13.39 cm)	29.5 in (11.61 cm)	29.5 in (11.61 cm)	34 in (13.39 cm)	34 in (13.39 cm)
SEAT BACK HEIGHT	29 in (73.7 cm)	28.75 in (73 cm)	29 in (73.7 cm)	28.75 in (73 cm)	29 in (73.7 cm)	28.75 in (73 cm)	29 in (73.7 cm)	28.75 in (73 cm)	29 in (73.7 cm)	28.75 in (73 cm)	29 in (73.7 cm)	28.75 in (73 cm)
OVERALL LENGTH	Upright: 39.5 in (100.3 cm) Reclined: 69 in (175.3 cm)	Upright: 41 in (104.1 cm) Reclined: 69.25 in (175.9 cm)	Upright: 39.5 in (100.3 cm) Reclined: 69 in (175.3 cm)	Upright: 41 in (104.1 cm) Reclined: 69.25 in (175.9 cm)	Upright: 39.5 in (100.3 cm) Reclined: 69 in (175.3 cm)	Upright: 41 in (104.1 cm) Reclined: 69.25 in (175.9 cm)	Upright: 39.5 in (100.3 cm) Reclined: 69 in (175.3 cm)	Upright: 41 in (104.1 cm) Reclined: 69.25 in (175.9 cm)	Upright: 39.5 in (100.3 cm) Reclined: 69 in (175.3 cm)	Upright: 41 in (104.1 cm) Reclined: 69.25 in (175.9 cm)	Upright: 39.5 in (100.3 cm) Reclined: 69 in (175.3 cm)	Upright: 41 in (104.1 cm) Reclined: 69.25 in (175.9 cm)
ARM HEIGHT	25.5 in (64.8 cm)	25.5 in (64.8 cm)	25.5 in (64.8 cm)	25.5 in (64.8 cm)	26" (66 cm)	26 in (66 cm)	26" (66 cm)	26 in (66 cm)	29 in (73.7 cm)	29 in (73.7 cm)	29 in (73.7 cm)	29 in (73.7 cm)
WEIGHT (NO OPTIONS)*	Push Back: 104 lb (47.17 kg) Gas Back: 108 lb (48.99 kg)	140 lb (63.5 kg)	Push Back: 122 lb (55.31 kg) Gas Back: 126 lb (57.15 kg)	158 lb (71.67 kg)	Push Back: 104 lb (47.17 kg) Gas Back: 108 lb (48.99 kg)	140 lb (63.5 kg)	Push Back: 122 lb (55.31 kg) Gas Back: 126 lb (57.15 kg)	158 lb (71.67 kg)	Push Back: 104 lb (47.17 kg) Gas Back: 108 lb (48.99 kg)	140 lb (63.5 kg)	Push Back: 122 lb (55.31 kg) Gas Back: 126 lb (57.15 kg)	158 lb (71.67 kg)
WEIGHT CAPACITY	350 lb (158.8 kg)	350 lb (158.8 kg)	500 lb (226.8 kg)	500 lb (226.8 kg)	350 lb (158.8 kg)	350 lb (158.8 kg)	500 lb (226.8 kg)	500 lb (226.8 kg)	350 lb (158.8 kg)	350 lb (158.8 kg)	500 lb (226.8 kg)	500 lb (226.8 kg)

## OPTIONS

	BATTERY BACKUP	DUAL SWITCHES	WIRED PENDANT	SWING ARMS	STANDARD TABLE <i>L, R or Both Sides</i>	PREMIUM TABLE <i>L, R or Both Sides</i>	FOOTPLATE	HEAT	HEAT & MASSAGE	IV POLE <i>L, R or Both Sides</i>	PUSH HANDLE
AVAILABLE FOR	POWER MODELS ONLY	POWER MODELS ONLY	POWER MODELS ONLY	ALL MODELS	ALL MODELS	ALL MODELS	5 in CASTER MODELS ONLY	ALL MODELS	ALL MODELS	ALL MODELS	<i>Recommended for Footplate models only</i>
*ADDITIONAL WEIGHT	-4.75 lb	-	-	-	+5 lb	+5 lb	-10 lb	-4 lb	-4 lb	+5 lb	-2 lb
*ADDITIONAL OVERALL WIDTH	-	-	-	Open 180°: -28.5 in (standard) / -32.5 in (XL)	Table down: -2.75 in (standard) / -2.5 in (XL) Table up: -13 in (standard) / -12.75 in (XL) <i>Measurements are per table</i>	Table down: -1.25 in (standard) / -1 in (XL) Table up: -11 in (standard) / -10.75 in (XL) <i>Measurements are per table</i>	-	-	-	-	-

Additional options are available. Please contact us to ask about our full line of options.

**WARRANTY** Contact Winco or visit our website for warranty details.

**COLORS** Contact Winco or visit our website for available color selections.

**NOTE** Winco reserves the right to change specifications without notice.

## CERTIFICATIONS

- 1) Winco recliners meet ANSI/BIFMA X5.4 for reliability testing
- 2) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture
- 3) NR or TB options are available. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for Upholstered Furniture
- 4) PRC models and Heat and Massage, Heat Only Options: Meet electrical safety standards, UL 1647 and CAN/ CSA-C22.2 No.68-92.
- 5) All individual electrical components for Heat and Massage / Heat Only options and power recline certified to IEC 60601-1.



## CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:  
**800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com**



## VERŌ PRC

POWER RECLINING CHAIR



### KEY FEATURES

- Sealed two-button power switch embedded on the side of the chair (dual switches and/or wired pendant options available)
- Easy to clean and service
- Durable construction - thoroughly tested to withstand the demands of healthcare environments
- Power-return from Trendelenburg (Trendelenburg models only)
- **4-year limited warranty** on the power recline control module and actuator motor

## Switch to the new Verō PRC

The Verō Power Reclining Chair (PRC) is the newest addition to our popular line of infinitely configurable Verō recliners. Verō PRC is comfortable and convenient for the patient while offering exceptional benefits to the caregiver and added value for facility managers. Caregivers can now assist patients up to 500 lbs at the push of a button with no fear of strain. Trendelenburg models are thoughtfully designed with the capability to **return patients upright from the trendelenburg position**, minimizing potential for caregiver injury.

A sealed power switch is embedded on the patient right side of the recliner for easy access. Add an optional second switch, at patient left side and/or a wired hand pendant for additional patient and caregiver convenience. Verō PRC powers the entire line of Verō recliners offering unmatched configurations, durability, patient comfort and caregiver peace of mind.

Contact Winco today to arrange for a demonstration of the Verō PRC or to request a quote!

## SPECIFICATIONS: VERÒ PRC

FOOT OPTION RECLINER SIZE	PEDESTAL FEET		3" DUAL CASTERS		5" CASTERS	
	STANDARD	X-LARGE	STANDARD	X-LARGE	STANDARD	X-LARGE
SEAT HEIGHT TO FLOOR	20.25 in (51.4 cm)		20.75 in (52.7 cm)		23.75 in (60.3 cm)	
SEAT WIDTH	22.5 in (57.2 cm)	26.5 in (67.3 cm)	22.5 in (57.2 cm)	26.5 in (67.3 cm)	22.5 in (57.2 cm)	26.5 in (67.3 cm)
SEAT DEPTH	23.5 in (59.7 cm)					
OVERALL HEIGHT	44 in (111.8 cm)		44.5 in (113 cm)		47.5 in (120.7 cm)	
OVERALL WIDTH (FIXED ARMS / NO TABLES)*	29.5 in (11.6 cm)	34 in (13.4 cm)	29.5 in (11.6 cm)	34 in (13.4 cm)	29.5 in (11.6 cm)	34 in (13.4 cm)
SEAT BACK HEIGHT	28.75 in ( 73 cm)					
OVERALL LENGTH	Upright: 41 in (104.1 cm) Reclined: 69.25 in (175.9 cm)					
ARM HEIGHT	25.5 in (64.8 cm)		26 in (66 cm)		29 in ( 73.7 cm)	
WEIGHT (NO OPTIONS)	140 lb (63.5 kg)	158 lb (71.7 kg)	140 lb (63.5 kg)	158 lb (71.7 kg)	140 lb (63.5 kg)	158 lb (71.7 kg)
WEIGHT CAPACITY	350 lb (158.8 kg)	500 lb (226.8 kg)	350 lb (158.8 kg)	500 lb (226.8 kg)	350 lb (158.8 kg)	500 lb (226.8 kg)

## OPTIONS: VERÒ PRC

- Wired Hand Pendant
- Dual 180 Degree Swing Arms
- Trendelenburg Positioning
- Standard or XL Seat Widths
- 3" Dual Casters, 5" Casters or Pedestal Feet
- Heat Only or Heat and Massage
- Packages Sliding, Weight Bearing Footplate
- Side Tables and/or Pivot Table
- Battery Back-Up for Recline Functions
- IV Pole
- Push Handle
- Accessory Hook
- Custom Armrest and Headrest Covers
- Urethane Arm Caps



+ Many more clinically useful options and accessories - Visit our website for more information

**WARRANTY** Contact Winco or visit our website for warranty details.

**COLORS** Contact Winco or visit our website for available color selections.

**NOTE** Winco reserves the right to change specifications without notice.

### CERTIFICATIONS

- 1) Winco recliners meet ANSI/BIFMA X5.4 for reliability testing
- 2) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture
- 3) NR or TB options are available. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for Upholstered Furniture
- 4) PRC models and Heat and Massage, Heat Only Options: Meet electrical safety standards, UL 1647 and CAN/ CSA-C22.2 No.68-92.
- 5) All individual electrical components for Heat and Massage / Heat Only options and power recline certified to IEC 60601-1.



## CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:

**800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com**

# AUGUSTINE

TREATMENT RECLINER

PREMIUM COLLECTION



Developed in conjunction with Oncology clinics, the Augustine was designed to be 100% patient-centric; to be comfortable and help put the patient at ease.

- Enhanced range of motion for optimum, patient selectable, positioning
- Spacious seat provides personal space and luxurious comfort
- Thoughtful features to ease chair movement and cleaning
- Unmatched stability and refined styling with concealed casters

Augustine demonstrates Winco's commitment to developing products that meet the needs of the people that live in and touch the recliners on a daily basis.

**ENHANCING THE PATIENT  
AND CAREGIVER EXPERIENCE.**



## SERENITY: PERSONAL SPACE & INDEPENDENCE

Often the medical treatment environment unavoidably impinges upon the patient's personal space. Augustine addresses this with:

- A low seat height to ease entry and exit
- A spacious 26" seat width for extra room to stretch out
- Easy access recline controls that allow patients to obtain desired positioning without assistance

When asked about their needs, patients were clear they wanted a chair that provided them with independence and a relaxing experience.

AUGUSTINE IS DESIGNED TO MIMIC FURNITURE THAT YOU WOULD BE PROUD TO HAVE IN YOUR OWN HOME



*The Central Front Braking option*

## CLINICIAN INPUT: DESIGN ENHANCEMENT

Facility managers requested an easier way to thoroughly clean the recliner. Once removed, the Quick-Release Seat provides full access to the inner sides of the recliner, backrest and seat for fast and thorough cleaning.

Simply pull-up at the front and slide forward. The seat is securely tethered so it won't get separated from the recliner. To replace simply line up the tabs, slide into place and push down – done!

QUICK RELEASE SEAT: NO NEED FOR ANY TOOLS AND RISK OF LOST PARTS



## CLINICAL UTILITY: PRACTICAL INNOVATION

Augustine has the clinical features required in a high use medical environment – casters for ease of movement, rigorously tested for durability and thoughtful features for ease of cleaning. These chairs bring dignity, comfort and grace back to the clinical setting without forfeiting the necessary clinical utility.

- Standard rear central caster lock for impressive stability
- Quick Release Seat is easily removed for efficient and thorough cleaning
- An optional Central Front Brake allows secure braking from the front of the chair, with no obtrusive levers or mechanisms to catch on or impede staff, patients or visitors.





**ANATOMY OF THE AUGUSTINE TREATMENT RECLINER**





## SPECIFICATIONS

MODEL	Augustine Treatment Recliner (5001)
Arm Style	Fixed
Casters	Front: 3" Fixed Rear: 3" Swivel Central Locking
Weight Capacity	500 lbs (226.80 kg)
Overall Height	45.75" (116.2 cm)
Overall Width*	34.5" (87.63 cm)
Seat Height (from Floor)	19.75" (50.2 cm)
Seat Width (Between Arms)	26" (66.04 cm)
Seat Depth	20.75" (52.71 cm)
Seat Back Height	28.5" (72.39 cm)
Armrest Height	26.75" (67.9 cm)
Upright Length	40.0" (101.6 cm)
Recline Length	76.0" (193.00 cm)
Range of Seat Angle	75° (Upright)   12° (Reclined)
Weight	181 lbs (82.1 kg)

OVER 50  
VINYL OPTIONS  
AVAILABLE



\* Add 2.5" to overall width for each table in the down position.

\* Add 11" to overall width for each raised side table.

All specifications subject to change without notice, please contact our customer care team to verify any critical dimensions prior to ordering.

## SUGGESTED OPTIONS & ACCESSORIES



### IV Pole (IV)

Left, Right or Both Sides



### Accessory Hook (AH)

Left, Right or Both Sides



### Heat Only (HT)



### Heat & Massage (HM)



### Premium Side Table (PT)

Left, Right or Both Sides



### Pivot Table (PV)

Left, Right or Both Sides



### Central Front Brake (BE)

**WARRANTY** Contact Winco or visit our website for warranty details. **COLORS** Contact Winco or visit our website for all available color selections. **NOTE** Winco reserves the right to change specifications without notice. **CERTIFICATIONS** 1) The Augustine Treatment Recliner meets ANSI/BIFMA X5.4 for reliability testing. 2) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture. 3) NR or TB options are available. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for upholstered furniture. 4) Heat and Massage and Heat Only Options meet electrical safety standards, UL 1647 and CAN/CSA-C22.2 No.68-92.



## CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:  
**800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com**

# XL CONVALESCENT RECLINER

5291 (WITH DUAL SIDE-TABLES)

5 SERIES



## Superior Comfort and Support...

- Built-in headrest
- Removable side panels for ease of cleaning
- 3 Locking positions: upright, TV or recline
- Dual side tables
- Padded armrests with scuff guards
- Recessed handgrips provide additional assistance for ingress and egress


ENHANCING THE PATIENT & CAREGIVER EXPERIENCE



# SPECIFICATIONS

DROP-ARM CONVALESCENT RECLINER	5291
DESCRIPTIONS	DUAL SIDE TABLES INCLUDED
Seat Height	22.5" (54.2 cm)
Seat Width	25" (63.5 cm)
Seat Depth	20" (50.8 cm)
Overall Height	49" (124.5 cm)
Overall Width (side tables down)	36.25" (90.1 cm)
Range of Seat Angle	8° - 32°
Range of Back Angle	62° - 38°
Back Height above Seat	31.5" (80 cm)
Overall Length (Back Upright)	44" (111.8 cm)
Overall Length (Back Reclined)	62" (157.5 cm)
Weight	112 lbs (50.8 kg)
Weight Capacity	450 lbs (204.1 kg)

## SUGGESTED OPTIONS & ACCESSORIES



**IV POLE**  
Left, Right, or  
Dual

OVER 50  
DESIGNER  
COLORS  
AVAILABLE

*Additional options and accessories available.*

*Please visit our website to view a full selection of options or contact your local representative.*

**WARRANTY** Contact Winco or visit our website for warranty details.

**COLORS** Visit our website or contact Winco for available color selections and/or swatches

**NOTE** Winco reserves the right to change specifications without notice.

**CERTIFICATIONS** 1) Winco recliners meet ANSI/BIFMA X5.4 for reliability testing 2) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture 3) NR or TB options are available. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for Upholstered Furniture 4) Heat and Massage, Heat Only Options: Meets electrical safety standards, UL 1647 and CAN/CSA-C22.2 No.68-92



## CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:

**800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com**

# PREMIER LIFECARE CLINER

5400

5 SERIES



## KEY FEATURES

- 3" Casters
- Urethane armrest covers
- Fixed arms
- 275lb Weight capacity
- 3-Position push-back recline
- Built-in patient positioning lock-bar

## OUTSTANDING VALUE AND QUALITY

Lower price does not mean lower quality when you purchase a Premier LifCare Cliner.

ENHANCING THE PATIENT & CAREGIVER EXPERIENCE



# ANATOMY OF A PREMIER LIFECARE CLINER

Full-coverage elastic support system in seat and frame back provide total support and enhanced comfort.

Patient-positioning lock-bar.

3" nylon casters for easy maneuverability.

Optional blow-molded, non-porous side tables. Optional wood upgrade also available.

Contoured urethane armrests  
\*wood options also available.



## SPECIFICATIONS

MODEL	5400
Seat Height	19.75"
Seat Width (Between Arms)	19.5"
Seat Depth	20"
Overall Height	45.25"
Overall Width*	27.75"
Back Height (Above Seat)	30.5"
Overall Length (Upright)	39.5"
Overall Length (Reclined)	62"
Arm Height	27.5"
Weight	90 lbs
Weight Capacity	275 lbs

## SUGGESTED OPTIONS & ACCESSORIES



**Accessory Hook (AH)**  
Left, Right or Both Sides



**Oxygen Tank Holder O2)**  
Left, Right or Both Sides



**Contoured Headrest (CHOO)**

\* Add 1" to chair width for each side table in down position. Add 10.5" to chair width for each raised side table.

**WARRANTY** Contact Winco or visit our website for warranty details.

**COLORS** Visit our website or contact Winco for available color selections and/or swatches

**NOTE** Winco reserves the right to change specifications without notice.

**CERTIFICATIONS** 1) Winco recliners meet ANSI/BIFMA X5.4 for reliability testing 2) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture 3) NR or TB options are available. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for Upholstered Furniture



## CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:

**800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com**

193915 REV A

# LIFE CARE CLINER

5 SERIES

5851 (WITH TRAY) / 5861 (WITHOUT TRAY)



## Value, Durability and Comfort...

- Molded side-panels
- Push handle for easy maneuverability
- 3 Locking positions: upright, TV or recline
- Padded armrests with scuff guards
- Recessed handgrips provide additional assistance for ingress and egress
- Model 5851: Side storage for lap tray

ENHANCING THE PATIENT & CAREGIVER EXPERIENCE



# SPECIFICATIONS

LIFE CARE CLINER	5851	5861
DESCRIPTIONS	TRAY INCLUDED	NO TRAY
Seat Height	20" (50.8 cm)	
Seat Width	19" (48.3 cm)	
Seat Depth	20" (50.8 cm)	
Overall Height	46.5" (118.1 cm)	
Overall Width	26.5" (67.3 cm)	
Range of Seat Angle	8° - 32°	
Range of Back Angle	62° - 38°	
Back Height above Seat	32" (81.3 cm)	
Overall Length (Back Upright)	41" (104.1 cm)	
Overall Length (Back Reclined)	63" (160 cm)	
Weight	85 lbs (38.6 kg)	80 lbs (36.3 kg)
Weight Capacity	275lbs (124.7 kg)	

## SUGGESTED OPTIONS & ACCESSORIES



**IV POLE**  
Left, Right, or  
Dual

OVER 50  
DESIGNER  
COLORS  
AVAILABLE

*Additional options and accessories available.*

*Please visit our website to view a full selection of options or contact your local representative.*

**WARRANTY** Contact Winco or visit our website for warranty details.

**COLORS** Visit our website or contact Winco for available color selections and/or swatches

**NOTE** Winco reserves the right to change specifications without notice.

**CERTIFICATIONS** 1) Winco recliners meet ANSI/BIFMA X5.4 for reliability testing 2) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture 3) NR or TB options are available. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for Upholstered Furniture 4) Heat and Massage, Heat Only Options: Meets electrical safety standards, UL 1647 and CAN/CSA-C22.2 No.68-92



## CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:

**800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com**





## TRENDELENBURG RECLINERS

6000 Series



# VERO

## TREATMENT RECLINER

PREMIUM SERIES



*Verō PRC with optional 3" dual casters, urethane arm caps and headrest cover & pillow set – shown in Bottle (85) and Limestone (57).*

Continuously evolving patient care requirements create opportunities for enhanced solutions. The Verō provides total flexibility in matching your application with exactly the features you require and the comfort your patients demand.

The Verō recliner's modular design allows clinicians to choose the best configuration to perfectly suit their specific needs:

- New Power Recline model (Verō PRC)!
- Multiple recline options and infinite position backrest promote patient independence.
- Bilateral and foot-operated controls provide clinician comfort and convenience.
- Durable construction and easy-to-clean surfaces enable increased efficiency.

**DESIGNED BY WINCO, BUILT BY YOU.  
THAT'S JUST SMART.**





# VERO CARE CLINER

PREMIUM SERIES

## ANATOMY OF THE VERO CARE CLINER

**Fastener-Free Headrest Cover**

*Optional for 5 Series*

**Push Handle**

*Option*

**Infinite Position Backrest**

**Power Recline**

*PRC Only*

**Gas Assist Control Lever**

*5X and 6X models only*

**Accessory Hook**

*Option*

**Power Recline Pendant**

*Option: PRC Only*

**Urethane Arm Cap**

*Option*

**Intuitive Side Tables**

*Option*

*Premium version shown*

**Foot Operated Swing Arm Release**

*Swing Arm Models only*

**Weight-Bearing, Retractable Foot Plate Option**

*5" Caster models only*

**Bilateral Trendelenburg Access**

*6 Series Models only*

**IV Pole + Holder**

*Option*

**Easily Accessed Heat & Massage or Heat Only Switch**

*Option*

50+ VINYL COLOR OPTIONS AVAILABLE

**Footrest with "Heel Drop" Comfort Zone**

**Independent Leg Rest**

**Durable, Easy-to-Clean, Powder-Coated Frame**

**180 Degree Swing Arm**

*Swing Arm models only*

**Fixed Front Casters with Accessible Lock Tabs**

*5" Casters only  
(Pedestal Feet or 3" Caster options available)*



## TRENDELENBURG

**Verō** Recliners are available with Trendelenburg (6 series) or without (5 series). On the **Verō**, Trendelenburg is controlled with bilateral access foot levers, reducing clinician's need to bend to reach controls. PRC models (with power recline) offer push-button recovery from Trendelenburg to minimize potential for clinician back strain.

Exclusive to **Verō**, Secure-Glide™ provides gentle patient positioning while avoiding injury to caregivers. The controlled descent helps eliminate any need or natural reflex to manually control the rapid patient descent in emergency situations.



Secure-Glide™ for stress-free positioning



Foot lever for exclusive Secure-Glide™ Trendelenburg descent


## RECLINE CONTROL OPTIONS

All **Verō** models feature convenient and comfortable infinite positioning of the back rest. With **Verō**, you have three options available for controlling recliner positioning:

**Power** Easy recline and recovery from trendelenburg at the push of a button. One switch standard. Dual switches and/or pendant options available.

**Gas Assist** provides positive locking position of the backrest angle with gas cylinder assisted operation. Locking action is controlled by easily accessed bilateral control levers.

**Push-Back** allows seamless positioning of the desired back rest angle. The gravity operated, gentle resistance designed into the chair maintains selected positioning.



**Power Recline**  
PRC models



**Gas Assist**  
5X and 6X models



**Push Back**  
5X and 6X models

## ARM OPTIONS

Functionality is enhanced through the availability of two arm configurations. Standard fixed arms are cost-effective, sturdy and durable. Optional swing arms allow 180 degree movement of the arm for convenient cleaning access and barrier-free patient transfer.

**Swing Arm Option** Swing arms are ideal when patient transfer and easy access are important considerations. Each side of the **Verō** is easily opened with convenient foot controls.

**Fixed Arm Option** Fixed arms are suitable for most general purpose operations when transfer of patient is not required or where working space is limited.



Upholstered Swing Arm Shown

## ARM CAP OPTIONS

**Urethane Cap Option** Urethane arm caps provides a unique, upscale appearance and added durability.

**Urethane Arm Cap**  
Option





## MATCH YOUR PURPOSE

The **Verō** modular design allows choices in how the chair is anchored or moved. Whether stationary or for patient transport\*, **Verō** has a perfect solution.

**Pedestal Feet** When mobility is not required or desired, pedestal feet provide a less clinical look and a lower seat height.

**Casters** provide ease of movement for cleaning, maintenance and staging. Two size options are available. Both include lockable, foot lever release controls.

**3" Casters** Lower seated height with the flexibility of mobility.

**5" Casters** are the optimal choice for ease in movement.



*It is recommended that after moving the chair always lock the casters prior to patient ingress/egress.*

## FOOT PLATE

Patients will appreciate the added comfort of the retractable foot plate.

The optional footplate is available for models with 5" casters. Its sturdy design is weight-bearing.

The footplate locks as weight is applied to prevent unwanted or unexpected movement.

**Note: Patients should be assisted at all times by a staff member when using the footplate.**



*Footplate is easily positioned by foot.*

## PUSH HANDLE OPTION

The push handle option is recommended in combination with footplate models if the chair is to be used for in-facility patient transport\*.

*\*Only 5X, 6X or PRC models with footplate and push handle are suitable for patient transport. See User Manual for important details.*



## SIZE / CAPACITY OPTIONS FOR MAXIMIZED FLEXIBILITY

The **Verō** is available in two sizes: standard and X-large. No matter which size is specified, **Verō** offers comfort, beauty and durability.

- Standard seat width is 22.5" between arms with a 350 pound weight capacity. This chair makes petite to average size patients feel right at home.
- X-large seat width is 26.5" between arms with a 500 pound weight capacity, making it suitable to care for Bariatric patients and provide a comfortable experience.



X-Large 500 lb Capacity



Standard 350 lb Capacity

## TABLE OPTIONS

There are three table options; standard, pivot and premium. All tables on the **Verō** are available in left, right or both configurations and can be ordered in a variety of combinations to suit facility requirements and clinician preference. All feature a stain resistant, non-porous material and all but the pivot table have a lip edge design to contain small spills.



Standard Table



Pivot Table




Premium Table

Each standard table adds 2.75" (2.5" for XL) to the overall chair width when lowered. Each premium table adds 1.25" (1" for XL) to the overall chair width when lowered. This added width should be considered when ordering table options appropriate for your facility.

- **Standard table** This value priced lift-to-raise-and-lower table will be familiar to current Winco Elite series owners.
- **Pivot table** swivels so patient can have the table in front of them.
- **Premium table** features an intuitive single hand up/down feature with audible positive lock and release.

## COMPLEMENT DESIGN & UTILITY

- **Heat Only Option** provides three patient-selected soothing warmth ranges for maximum comfort.
- **Heat and Massage Option** adds a selection of five gentle vibrating massage modes to the heat option, enhancing patient comfort level.
- **Optional Accessory Hook** The stainless steel hook remains flush until needed, making it a convenient choice for a Foley bag. Available for either left, right or both sides.
- **Optional IV Pole and Holder** for either left, right or both sides.
- **Headrest Cover and Armrest Covers** Covers extend the life of the headrest and armrests and are fastener-free, making them easy to lift on and off.



**Heat & Massage**  
3 levels of heat, 5 levels of massage

**Heat Only**  
Low/Med/High

CONVENIENTLY  
PLACED CONTROLS  
FOR PATIENT  
INDEPENDENCE



**Accessory Hook**  
Shown here with Foley bag



**IV Holder Mount**  
Available for left, right or both sides



# MODELS & SPECIFICATIONS

MOBILITY OPTION	PEDESTAL FEET				3" CASTERS				5" CASTERS			
	STANDARD		X-LARGE		STANDARD		X-LARGE		STANDARD		X-LARGE	
RECLINER SIZE												
MODEL OPTION	PUSH BACK/ GAS BACK	POWER	PUSH BACK/ GAS BACK	POWER	PUSH BACK/ GAS BACK	POWER	PUSH BACK/ GAS BACK	POWER	PUSH BACK/ GAS BACK	POWER	PUSH BACK/ GAS BACK	POWER
SEAT HEIGHT TO FLOOR	19.25 in (48.9 cm)	20.25 in (51.4 cm)	19.25 in (48.9 cm)	20.25 in (51.4 cm)	19.75 in (50.2 cm)	20.75 in (52.7 cm)	19.75 in (50.2 cm)	20.75 in (52.7 cm)	22.75 in (57.8 cm)	23.75 in (60.3 cm)	22.75 in (57.8 cm)	23.75 in (60.3 cm)
SEAT WIDTH	22.5 in (57.2 cm)	22.5 in (57.2 cm)	26.5 in (67.3 cm)	26.5 in (67.3 cm)	22.5 in (57.2 cm)	22.5 in (57.2 cm)	26.5 in (67.3 cm)	26.5 in (67.3 cm)	22.5 in (57.2 cm)	22.5 in (57.2 cm)	26.5 in (67.3 cm)	26.5 in (67.3 cm)
SEAT DEPTH	22 in (55.9 cm)	23.5 in (59.7 cm)	22 in (55.9 cm)	23.5 in (59.7 cm)	22 in (55.9 cm)	23.5 in (59.7 cm)	22 in (55.9 cm)	23.5 in (59.7 cm)	22 in (55.9 cm)	23.5 in (59.7 cm)	22 in (55.9 cm)	23.5 in (59.7 cm)
OVERALL HEIGHT	44.75 in (113.7 cm)	44 in (111.8 cm)	44.75 in (113.7 cm)	44 in (111.8 cm)	45.25 in (115 cm)	44.5 in (113 cm)	45.25 in (115 cm)	44.5 in (113 cm)	48.25 in (122.6 cm)	47.5 in (120.7 cm)	48.25 in (122.6 cm)	47.5 in (120.7 cm)
OVERALL WIDTH (FIXED ARMS / NO TABLES)*	29.5 in (11.61 cm)	29.5 in (11.61 cm)	34 in (13.39 cm)	34 in (13.39 cm)	29.5 in (11.61 cm)	29.5 in (11.61 cm)	34 in (13.39 cm)	34 in (13.39 cm)	29.5 in (11.61 cm)	29.5 in (11.61 cm)	34 in (13.39 cm)	34 in (13.39 cm)
SEAT BACK HEIGHT	29 in (73.7 cm)	28.75 in (73 cm)	29 in (73.7 cm)	28.75 in (73 cm)	29 in (73.7 cm)	28.75 in (73 cm)	29 in (73.7 cm)	28.75 in (73 cm)	29 in (73.7 cm)	28.75 in (73 cm)	29 in (73.7 cm)	28.75 in (73 cm)
OVERALL LENGTH	Upright: 39.5 in (100.3 cm) Reclined: 69 in (175.3 cm)	Upright: 41 in (104.1 cm) Reclined: 69.25 in (175.9 cm)	Upright: 39.5 in (100.3 cm) Reclined: 69 in (175.3 cm)	Upright: 41 in (104.1 cm) Reclined: 69.25 in (175.9 cm)	Upright: 39.5 in (100.3 cm) Reclined: 69 in (175.3 cm)	Upright: 41 in (104.1 cm) Reclined: 69.25 in (175.9 cm)	Upright: 39.5 in (100.3 cm) Reclined: 69 in (175.3 cm)	Upright: 41 in (104.1 cm) Reclined: 69.25 in (175.9 cm)	Upright: 39.5 in (100.3 cm) Reclined: 69 in (175.3 cm)	Upright: 41 in (104.1 cm) Reclined: 69.25 in (175.9 cm)	Upright: 39.5 in (100.3 cm) Reclined: 69 in (175.3 cm)	Upright: 41 in (104.1 cm) Reclined: 69.25 in (175.9 cm)
ARM HEIGHT	25.5 in (64.8 cm)	25.5 in (64.8 cm)	25.5 in (64.8 cm)	25.5 in (64.8 cm)	26" (66 cm)	26 in (66 cm)	26" (66 cm)	26 in (66 cm)	29 in (73.7 cm)	29 in (73.7 cm)	29 in (73.7 cm)	29 in (73.7 cm)
WEIGHT (NO OPTIONS)*	Push Back: 104 lb (47.17 kg) Gas Back: 108 lb (48.99 kg)	140 lb (63.5 kg)	Push Back: 122 lb (55.31 kg) Gas Back: 126 lb (57.15 kg)	158 lb (71.67 kg)	Push Back: 104 lb (47.17 kg) Gas Back: 108 lb (48.99 kg)	140 lb (63.5 kg)	Push Back: 122 lb (55.31 kg) Gas Back: 126 lb (57.15 kg)	158 lb (71.67 kg)	Push Back: 104 lb (47.17 kg) Gas Back: 108 lb (48.99 kg)	140 lb (63.5 kg)	Push Back: 122 lb (55.31 kg) Gas Back: 126 lb (57.15 kg)	158 lb (71.67 kg)
WEIGHT CAPACITY	350 lb (158.8 kg)	350 lb (158.8 kg)	500 lb (226.8 kg)	500 lb (226.8 kg)	350 lb (158.8 kg)	350 lb (158.8 kg)	500 lb (226.8 kg)	500 lb (226.8 kg)	350 lb (158.8 kg)	350 lb (158.8 kg)	500 lb (226.8 kg)	500 lb (226.8 kg)

## OPTIONS

	BATTERY BACKUP	DUAL SWITCHES	WIRED PENDANT	SWING ARMS	STANDARD TABLE <i>L, R or Both Sides</i>	PREMIUM TABLE <i>L, R or Both Sides</i>	FOOTPLATE	HEAT	HEAT & MASSAGE	IV POLE <i>L, R or Both Sides</i>	PUSH HANDLE
AVAILABLE FOR	POWER MODELS ONLY	POWER MODELS ONLY	POWER MODELS ONLY	ALL MODELS	ALL MODELS	ALL MODELS	5 in CASTER MODELS ONLY	ALL MODELS	ALL MODELS	ALL MODELS	<i>Recommended for Footplate models only</i>
*ADDITIONAL WEIGHT	-4.75 lb	-	-	-	+5 lb	+5 lb	-10 lb	-4 lb	-4 lb	+5 lb	-2 lb
*ADDITIONAL OVERALL WIDTH	-	-	-	Open 180°: -28.5 in (standard) / -32.5 in (XL)	Table down: -2.75 in (standard) / -2.5 in (XL) Table up: -13 in (standard) / -12.75 in (XL) <i>Measurements are per table</i>	Table down: -1.25 in (standard) / -1 in (XL) Table up: -11 in (standard) / -10.75 in (XL) <i>Measurements are per table</i>	-	-	-	-	-

Additional options are available. Please contact us to ask about our full line of options.

**WARRANTY** Contact Winco or visit our website for warranty details.

**COLORS** Contact Winco or visit our website for available color selections.

**NOTE** Winco reserves the right to change specifications without notice.

## CERTIFICATIONS

- 1) Winco recliners meet ANSI/BIFMA X5.4 for reliability testing
- 2) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture
- 3) NR or TB options are available. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for Upholstered Furniture
- 4) PRC models and Heat and Massage, Heat Only Options: Meet electrical safety standards, UL 1647 and CAN/ CSA-C22.2 No.68-92.
- 5) All individual electrical components for Heat and Massage / Heat Only options and power recline certified to IEC 60601-1.



## CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:  
**800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com**



## VERŌ PRC

POWER RECLINING CHAIR



### KEY FEATURES

- Sealed two-button power switch embedded on the side of the chair (dual switches and/or wired pendant options available)
- Easy to clean and service
- Durable construction - thoroughly tested to withstand the demands of healthcare environments
- Power-return from Trendelenburg (Trendelenburg models only)
- **4-year limited warranty** on the power recline control module and actuator motor

## Switch to the new Verō PRC

The Verō Power Reclining Chair (PRC) is the newest addition to our popular line of infinitely configurable Verō recliners. Verō PRC is comfortable and convenient for the patient while offering exceptional benefits to the caregiver and added value for facility managers. Caregivers can now assist patients up to 500 lbs at the push of a button with no fear of strain. Trendelenburg models are thoughtfully designed with the capability to **return patients upright from the trendelenburg position**, minimizing potential for caregiver injury.

A sealed power switch is embedded on the patient right side of the recliner for easy access. Add an optional second switch, at patient left side and/or a wired hand pendant for additional patient and caregiver convenience. Verō PRC powers the entire line of Verō recliners offering unmatched configurations, durability, patient comfort and caregiver peace of mind.

Contact Winco today to arrange for a demonstration of the Verō PRC or to request a quote!

## SPECIFICATIONS: VERÒ PRC

FOOT OPTION RECLINER SIZE	PEDESTAL FEET		3" DUAL CASTERS		5" CASTERS	
	STANDARD	X-LARGE	STANDARD	X-LARGE	STANDARD	X-LARGE
SEAT HEIGHT TO FLOOR	20.25 in (51.4 cm)		20.75 in (52.7 cm)		23.75 in (60.3 cm)	
SEAT WIDTH	22.5 in (57.2 cm)	26.5 in (67.3 cm)	22.5 in (57.2 cm)	26.5 in (67.3 cm)	22.5 in (57.2 cm)	26.5 in (67.3 cm)
SEAT DEPTH	23.5 in (59.7 cm)					
OVERALL HEIGHT	44 in (111.8 cm)		44.5 in (113 cm)		47.5 in (120.7 cm)	
OVERALL WIDTH (FIXED ARMS / NO TABLES)*	29.5 in (11.6 cm)	34 in (13.4 cm)	29.5 in (11.6 cm)	34 in (13.4 cm)	29.5 in (11.6 cm)	34 in (13.4 cm)
SEAT BACK HEIGHT	28.75 in ( 73 cm)					
OVERALL LENGTH	Upright: 41 in (104.1 cm) Reclined: 69.25 in (175.9 cm)					
ARM HEIGHT	25.5 in (64.8 cm)		26 in (66 cm)		29 in ( 73.7 cm)	
WEIGHT (NO OPTIONS)	140 lb (63.5 kg)	158 lb (71.7 kg)	140 lb (63.5 kg)	158 lb (71.7 kg)	140 lb (63.5 kg)	158 lb (71.7 kg)
WEIGHT CAPACITY	350 lb (158.8 kg)	500 lb (226.8 kg)	350 lb (158.8 kg)	500 lb (226.8 kg)	350 lb (158.8 kg)	500 lb (226.8 kg)

## OPTIONS: VERÒ PRC

- Wired Hand Pendant 
- Dual 180 Degree Swing Arms
- Trendelenburg Positioning
- Standard or XL Seat Widths
- 3" Dual Casters, 5" Casters or Pedestal Feet
- Heat Only or Heat and Massage
- Packages Sliding, Weight Bearing Footplate
- Side Tables and/or Pivot Table
- Battery Back-Up for Recline Functions
- IV Pole
- Push Handle
- Accessory Hook
- Custom Armrest and Headrest Covers
- Urethane Arm Caps

+ Many more clinically useful options and accessories - Visit our website for more information

**WARRANTY** Contact Winco or visit our website for warranty details.

**COLORS** Contact Winco or visit our website for available color selections.

**NOTE** Winco reserves the right to change specifications without notice.

### CERTIFICATIONS

- 1) Winco recliners meet ANSI/BIFMA X5.4 for reliability testing
- 2) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture
- 3) NR or TB options are available. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for Upholstered Furniture
- 4) PRC models and Heat and Massage, Heat Only Options: Meet electrical safety standards, UL 1647 and CAN/ CSA-C22.2 No.68-92.
- 5) All individual electrical components for Heat and Massage / Heat Only options and power recline certified to IEC 60601-1.



## CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:

**800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com**

# INVERNESS

## 24-HOUR TREATMENT RECLINER

PREMIUM SERIES



The Inverness 24-Hour Treatment Recliner™ is superbly crafted with a suite of intuitive, practical and durable innovations to address patient and clinician comfort.

- Promotes patient independence with ultra-low effort, bilateral controls
- Supports clinician comfort with front-locking brake tabs and back-relieving side table heights
- Enhances facility efficiency with attractive, durable construction that is easy to clean

Thousands of documented hours of feedback and refinements create an ideal blend of daytime luxury with nocturnal comfort upon which patients truly can sleep through the night. The Inverness provides comfort, utility and flexibility for multiple treatment applications.

**ENHANCING THE PATIENT  
AND CAREGIVER EXPERIENCE.**



## DESIGNED WITH PATIENTS

Patient independence, with no need to rely on staff for assistance. We listened to patients and included:

- An infinitely adjustable backrest that locks in the chosen position
- An innovative design requiring 50% less pull force for one-handed positioning
- Patient comfort in every position; full upright to reclined
- Wide padded seats, naturally curved armrests, accessible controls, and strategically located lumbar support

During development, many patients noted that they would choose clinics with this chair over those with another recliner model.

“I REALLY LOVED THE CHAIR AND WAS ABLE TO RELAX AND GET TREATMENT IN COMFORT.” - *William, Patient*



## DESIGNED WITH THE FACILITY

The Inverness 24-Hour Treatment Recliner's form and function make it a practical and valuable tool, easily expanding nocturnal capacity to existing clinics.

- Both arms swing out 180 degrees for barrier free patient transfer
- Innovative safety shields direct errant fluids and make cleaning a snap.
- Generous seat width and 500 pound weight capacity eliminate the need for multiple models.
- Swing arm design allows deployment of Trendelenburg function even when fully open.

We also took into account the unique requirements of older retrofit clinics. With the arms fully open, the Inverness will easily roll through a 34.5" opening!

“STAYS LOCKED SECURELY AND EASILY. CURVED ARM RESTS MORE COMFORTABLE FOR PATIENTS. LOOKS LESS CLINICAL. EASY TO USE AND MOVE. LESS PATIENT COMPLAINTS. PATIENTS ARE ABLE TO EASILY ADJUST BACK DURING TREATMENT.” - *Bridgette D., Clinic Manager*



## DESIGNED WITH CLINICANS

We also listened to the clinician and designed:

- Higher patient positioning with side tables at an ideal level to reduce potential for back strain.
- Bilateral controls that are easily accessed from either side.
- Foot operated swing-arm latch releases, unique foot plate lock and Trendelenburg controls to reduce bending, stooping and straining.
- Gas cylinder and spring supported mechanism to significantly ease any required effort to help patient recline or sit upright.

We've strived to address everything, including never hidden, always accessible, front caster locks.

“GREAT CHAIR, EASY ON TECH'S BACK. BETTER OVERALL FOR THE PATIENT. EASY TO CLEAN AND MOVE. NICE JOB, GUYS AT WINCO! THANK YOU FOR THINKING OF THE WORKERS ALSO!” - *Christopher F., Technician*





**ANATOMY OF THE INVERNESS 24-HOUR TREATMENT RECLINER**

**Fastener-Free Headrest Cover**

*Option*

**Push Handle**

**Infinite Position Backrest**

**Easy-Pull Bilateral Recline Handles**

**Accessory Hook**

*Option*

**Bilateral Trendelenburg Access**

**Footrest with "Heel Drop" Comfort Zone**

**Intuitive "Comfort-Height" Side Tables**

*Option*

**Foot-Operated Arm Release**

**IV Pole & Holder**

*Option*

**Easily Accessed Heat & Massage Switch**

*Option*

**Independent Leg Rest**

**Safety Shield Mechanism Guards**

**Durable, Easy-to-Clean Powder-Coated Frame**

**Locking, Retractable Foot Plate**

**Fixed Front Casters with Accessible Lock Tabs**

**180° Swing Arm**



## SPECIFICATIONS

Arm Style	Dual 180° Swing Arms
Casters	Front: 5" fixed, locking Rear: 5" Swivel locking
Weight Capacity	500 lbs (226.80 kg)
Overall Height	51.75" (131.45 cm)
Overall Width*	34.5" (87.63 cm)
Seat Height (from Floor)	25.5" (64.77 cm)
Seat Width (Between Arms)	26" (66.04 cm)
Seat Depth	20.75" (52.71 cm)
Seat Back Height	28.5" (72.39 cm)
Armrest Height	32.5" (82.55 cm)
Upright Length	41.5" (105.41 cm)
Recline Length	75.5" (191.77 cm)
Arm-Span (Swing Arms)	105.75" (268.61 cm)
Weight	199 lbs (90.26 kg)

OVER 50  
VINYL OPTIONS  
AVAILABLE



\* Add 2" to overall width for each table in the down position.

\* Add 11" to overall width for each raised side table.

All specifications subject to change without notice, please contact our customer care team to verify any critical dimensions prior to ordering.

## SUGGESTED OPTIONS & ACCESSORIES



**IV Pole (IV)**  
Left, Right  
or Both Sides



**Accessory Hook (AH)**  
Left, Right  
or Both Sides



**Heat Only (HT)**



**Heat & Massage (HM)**



**Premium Side Table (PT)**  
Left, Right  
or Both Sides



**Pivot Table (PV)**  
Left, Right  
or Both Sides

**WARRANTY** Contact Winco or visit our website for warranty details. **COLORS** Contact Winco or visit our website for all available color selections. **NOTE** Winco reserves the right to change specifications without notice. **CERTIFICATIONS** 1) The Inverness 24-Hour Treatment Recliner™ meets ANSI/BIFMA X5.4 for reliability testing 2) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture 3) NR or TB options are available. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for upholstered furniture. 4) Heat and Massage, Heat Only Options: Meets electrical safety standards, UL 1647 and CAN/CSA-C22.2 No.68-92.



## CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:  
**800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com**



# CARE CLINER

6530 (STANDARD, NYLON CASTERS) / 6531 (STANDARD, STEEL CASTERS)  
6540 (XL, NYLON CASTERS) / 6541 (XL, STEEL CASTERS)

6 SERIES



## KEY FEATURES

- Standard headrest cover
- Padded armrests
- Dual fold-down side tables
- Removable side panels
- Trendelenburg positioning

## OUTSTANDING VALUE AND QUALITY

Lower price does not mean lower quality when you purchase a Winco Care Cliner.

Our Care Cliner comes with all the superior features you might find on higher-priced clinical chairs

ENHANCING THE PATIENT & CAREGIVER EXPERIENCE



# ANATOMY OF A CARE CLINER

Trendelenberg release on back with easy-access handle .

Recessed hand grips to assist chair ingress and egress

Non-porous, fold-down side tray contains small spills. Standard on right and left side.

Total-lock and directional-lock casters provide ease of maneuverability. \*Available in nylon or steel.

Removable side panels for easy cleaning and maintenance



## SPECIFICATIONS

MODELS	6530 / 6531	6540 / 6541
DESCRIPTIONS	STANDARD SIZE	XL SIZE
Seat Height	21" (53.3 cm)	21" (53.3 cm)
Seat Width	21" (53.3 cm)	25" (63.5 cm)
Seat Depth	20" (51 cm)	20" (51 cm)
Overall Height	49" (124.5 cm)	49" (124.5 cm)
Overall Width*	32.25" (81.9 cm)	36.25" (90.1 cm)
Back Height above Seat	31.5" (80 cm)	31.5" (80 cm)
Overall Length (Upright)	39" (99 cm)	44" (111.8 cm)
Overall Length (Reclined)	76.5" (194.3 cm)	76.5" (194.3 cm)
Armrest Height	27" (68.6 cm)	27" (68.6 cm)
Weight	94 lbs (42.6 kg)	111 lbs (50.4 kg)
Weight Capacity	275 lbs (124.7 kg)	450 lbs (204 kg)

\* Add 9.25" to Overall Chair Width for each raised side table.


## OPTIONS



Leg Extensions (RE)



Side Cushions (SC)



IV Pole (IV)



Heat Only (H)



Heat and Massage (HM)

**WARRANTY** Contact Winco or visit our website for warranty details.

**COLORS** Visit our website or contact Winco for available color selections and/or swatches

**NOTE** Winco reserves the right to change specifications without notice.

**CERTIFICATIONS** 1) Winco recliners meet ANSI/BIFMA X5.4 for reliability testing 2) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture 3) NR or TB options are available. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for Upholstered Furniture 4) Heat and Massage, Heat Only Options: Meets electrical safety standards, UL 1647 and CAN/CSA-C22.2 No.68-92



## CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:

800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com

# DROP-ARM CARE CLINER

6550 (STANDARD, NYLON CASTERS) / 6551 (STANDARD, STEEL CASTERS)  
6570 (XL, NYLON CASTERS) / 6571 (XL, STEEL CASTERS)

6 SERIES



## KEY FEATURES

- Standard headrest cover
- Padded drop-down armrests
- Dual fold-down side tables
- Easy-clean side panels
- Trendelenburg positioning

6570 Shown

## VERSATILITY, VALUE AND QUALITY

Our Drop-Arm Care Cliner comes with all the superior features you might find on higher-priced clinical chairs and make it easier for caregivers to assist their patients!

ENHANCING THE PATIENT & CAREGIVER EXPERIENCE



# ANATOMY OF A DROP-ARM CARE CLINER

Trendelenberg release on back with easy-access handle .

Non-porous, fold-down side tray contains small spills. Standard on right and left side.

One-touch release drop-arms

Total-lock and directional-lock casters provide ease of maneuverability. \*Available in nylon or steel.



## SPECIFICATIONS

MODELS	6550 / 6551	6570 / 6571
DESCRIPTIONS	STANDARD SIZE	XL SIZE
Seat Height	21" (53.3 cm)	21" (53.3 cm)
Seat Width	21" (53.3 cm)	25" (63.5 cm)
Seat Depth	20" (51 cm)	20" (51 cm)
Overall Height	49" (124.5 cm)	49" (124.5 cm)
Overall Width*	33.75" (85.73 cm)	37.75" (95.9 cm)
Back Height above Seat	31.5" (80 cm)	31.5" (80 cm)
Overall Length (Upright)	40" (101.6 cm)	43" (109.2 cm)
Overall Length (Reclined)	76.5" (194.3 cm)	76.5" (194.3 cm)
Armrest Height	27" (68.6 cm)	27" (68.6 cm)
Weight	105 lbs (47 kg)	111 lbs (50.4 kg)
Weight Capacity	275 lbs (124.7 kg)	450 lbs (204 kg)

\* Add 9.25" to Overall Chair Width for each raised side table.


## OPTIONS



Leg Extensions (SE)



Side Cushions (SC)



IV Pole (IV)



Heat Only (H)



Heat and Massage (HM)

**WARRANTY** Contact Winco or visit our website for warranty details.

**COLORS** Visit our website or contact Winco for available color selections and/or swatches

**NOTE** Winco reserves the right to change specifications without notice.

**CERTIFICATIONS** 1) Winco recliners meet ANSI/BIFMA X5.4 for reliability testing 2) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture 3) NR or TB options are available. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for Upholstered Furniture 4) Heat and Massage, Heat Only Options: Meets electrical safety standards, UL 1647 and CAN/CSA-C22.2 No.68-92



## CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:

800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com

# ELITE CARE CLINER

6900 (STANDARD) / 6910 (X-LARGE)

6 SERIES



6900 Shown

## Comfort, Durability and Value...

- Dual fold-down side tables
- Built-in headrest with cover
- Trendelenburg positioning
- 5" Total-lock and directional-lock casters

ENHANCING THE PATIENT & CAREGIVER EXPERIENCE



# SPECIFICATIONS

ELITE CARE CLINER	6900	6910
DESCRIPTIONS	STANDARD	X-LARGE
Seat Height	21" (53.3 cm)	
Seat Width	21" (53.3 cm)	25" (63.5 cm)
Seat Depth	20" (51 cm)	
Overall Height	48.5" (123.2 cm)	
Overall Width*	33.5" (85.1 cm)	37.5" (95.3 cm)
Range of Seat Angle	8° - 32°	
Range of Back Angle	62° - 38°	
Back Height above Seat	32.5" (82.6 cm)	
Overall Length (Back Upright)	41" (111.8 cm)	
Overall Length (Back Reclined)	76.5" (194.3 cm)	
Weight	100 lbs (45.4 kg)	123 lbs (55.8 kg)
Weight Capacity	350 lbs (158.7 kg)	450 lbs (204 kg)

\* Add 9" to overall chair width for each raised side-table.

## SUGGESTED OPTIONS & ACCESSORIES



**HEAT**  
(H)



**HEAT & MASSAGE**  
(H/M)



**IV POLE**  
Left, Right, or  
Dual

OVER 50  
DESIGNER  
COLORS  
AVAILABLE

*Additional options and accessories available.*

*Please visit our website to view a full selection of options or contact your local representative.*

**WARRANTY** Contact Winco or visit our website for warranty details.

**COLORS** Visit our website or contact Winco for available color selections and/or swatches

**NOTE** Winco reserves the right to change specifications without notice.

**CERTIFICATIONS** 1) Winco recliners meet ANSI/BIFMA X5.4 for reliability testing 2) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture 3) NR or TB options are available. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for Upholstered Furniture 4) Heat and Massage, Heat Only Options: Meets electrical safety standards, UL 1647 and CAN/CSA-C22.2 No.68-92



## CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:

**800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com**



# SWING-ARM ELITE CARE CLINER

6 SERIES

6940 (STANDARD) / 6950 (X-LARGE)



## Ease Patient Transfers Without Sacrificing Patient Comfort...

- Dual swing-arms with foot-activated release levers
- Dual side-tables
- Built-in headrest with cover
- Trendelenburg positioning
- 5" Total-lock and directional-lock casters

ENHANCING THE PATIENT & CAREGIVER EXPERIENCE



# SPECIFICATIONS

SWING-ARM ELITE CARE CLINER	6940	6950
DESCRIPTIONS	STANDARD	X-LARGE
Seat Height	21" (53.3 cm)	
Seat Width	21" (53.3 cm)	25" (63.5 cm)
Seat Depth	20" (51 cm)	
Overall Height	48.5" (123.2 cm)	
Overall Width*	33.5" (85.1 cm)	37.5" (95.3 cm)
Range of Seat Angle	8° - 32°	
Range of Back Angle	62° - 38°	
Back Height above Seat	32.5" (82.6 cm)	
Overall Length (Back Upright)	41" (111.8 cm)	
Overall Length (Back Reclined)	76.5" (194.3 cm)	
Weight	104 lbs (47.2 kg)	127 lbs (57.6 kg)
Weight Capacity	350 lbs (158.7 kg)	450 lbs (204 kg)

\* Add 9" to overall chair width for each raised side-table.


## SUGGESTED OPTIONS & ACCESSORIES



**HEAT**  
(H)



**HEAT & MASSAGE**  
(H/M)



**IV POLE**  
Left, Right, or  
Dual

OVER 50  
DESIGNER  
COLORS  
AVAILABLE

*Additional options and accessories available.*

*Please visit our website to view a full selection of options or contact your local representative.*

**WARRANTY** Contact Winco or visit our website for warranty details.

**COLORS** Visit our website or contact Winco for available color selections and/or swatches

**NOTE** Winco reserves the right to change specifications without notice.

**CERTIFICATIONS** 1) Winco recliners meet ANSI/BIFMA X5.4 for reliability testing 2) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture 3) NR or TB options are available. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for Upholstered Furniture 4) Heat and Massage, Heat Only Options: Meets electrical safety standards, UL 1647 and CAN/CSA-C22.2 No.68-92



## CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:

**800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com**



## PHLEBOTOMY COLLECTION

Blood-Draw Chairs





*Shown: Solace BDC Model 2573-17*

BDC COLLECTION

ENHANCING THE PATIENT & CAREGIVER EXPERIENCE





# PHLEBOTOMY CHAIR HIGHLIGHTS

BDC COLLECTION

## FUNCTIONALITY, COMFORT & CONVENIENCE FOR PATIENTS AND STAFF

Multiple options to meet your application and patient treatment needs.

*\*Features and options vary by model*



*Shown: Solace BDC Model 2573-17*

ALL MODELS INCLUDE A LIFETIME FRAME WARRANTY FOR YEARS OF TROUBLE-FREE USE





## **SPIRIT BDC**

**2570 / 2571**

This basic model features two infinitely adjustable armrests, durable plastic resin seat and non-marring leveling feet. A clean and simple design that stands up to years of use.



## **SOLACE BDC**

**2572 / 2573**

Winco's Solace BDC features two infinitely adjustable armrests, an upholstered seat and back and non-marring leveling feet. The flip-up contoured urethane armrest provides comfort, easy access, and added security.



## **UNITY BDC**

**2574 / 2575**

This XL Phlebotomy chair features two infinitely adjustable armrests, an upholstered 24" seat and back and a 450lb. weight capacity.



## **HARMONY BDC**

**2580**

Featuring a 36.75" seat, an 800lb weight capacity and an array of options, the Harmony BDC can accommodate nearly any medical setting. Perfect for bariatric applications and pediatric facilities. Shown with dual flip arms. Additional literature available - visit our website for more details.

## AT-A-GLANCE MODEL GUIDE

MODEL		SEAT & BACK	ARMS					WEIGHT CAPACITY	FOOTREST
			FLIP-UP ARMREST	PIVOT ARM	SHORT ARMS	L-ARMS	CABINET		
SPIRIT BDC	2570	Resin	✓	●	N/A	N/A	✓	300 lbs.	N/A
	2571	Resin	✓	●	N/A	N/A	N/A	300 lbs.	N/A
SOLACE BDC	2572	Padded	✓	●	N/A	N/A	✓	300 lbs.	N/A
	2573	Padded	✓	●	N/A	N/A	N/A	300 lbs.	N/A
UNITY BDC	2574	Padded	✓	●	N/A	N/A	✓	450 lbs.	N/A
	2575	Padded	✓	●	N/A	N/A	N/A	450 lbs.	N/A
HARMONY BDC	2580	Padded	N/A	N/A	✓	●	●	800 lbs	●

✓ = Standard ● = Optional N/A= Not Available

## SUGGESTED FURNISHINGS TO COMPLETE YOUR EXAM ROOM



### Caregiver Task Seating

Deluxe Gas-Lift Task Stool  
Shown: Model 4400-17



### Privess™ Privacy Screens

Privess™ Elite Privacy Screen  
Shown: Model 3170-81



### Treatment Tables

Flat Top Treatment Table  
Shown: Model 8500-08

**WARRANTY** Contact Winco or visit our website for complete warranty details.

**VINYL COLORS** Contact Winco or visit our website for available color options.

**CERTIFICATIONS** 1) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture  
3) NR or TB options are available. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for upholstered furniture.



## CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:  
**800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com**

# HARMONY BDC

EXTRA WIDE BLOOD-DRAWING CHAIR

PHLEBOTOMY COLLECTION



*Dual L-Arm Option Shown*

Winco's Extra Wide Blood Drawing Chair combines functional design with quality and unmatched durability. The new Harmony BDC offers comfort and convenience with a variety of accessories to accommodate any medical setting where samples are drawn.

- Harmony BDC offered with multiple options to meet every clinic need.
- Model features include an adjustable-height, padded L-arm option for easy access and exceptional patient and clinician comfort.
- Industry leading 800 lb. capacity suitable for most patients or parent/child combinations.
- Six (6) high-impact, adjustable leveling feet provide stable seating on practically any surface.
- Manufactured with a sturdy, powder-coated steel frame for maximum strength and years of trouble-free service backed with a lifetime frame warranty.

**ENHANCING THE PATIENT & CAREGIVER EXPERIENCE**



## Designed for bariatric applications, yet perfect for pediatric settings.


In Family Practice or Pediatric Care, the Harmony BDC doubles as a one-size-fits-all patients chair or to safely and comfortably seat both parent and child, helping to soothe patient anxiety. The Harmony BDC represents one of the widest and highest weight capacity blood-draw chairs in the industry, tested and verified to BIFMA\* standards.



*Dual L-Arm and Footrest Options Shown*

### WINCO HAS A HISTORY OF DESIGNING AND MANUFACTURING ROBUST, THOUGHTFUL PHLEBOTOMY SEATING SOLUTIONS. HARMONY BDC FEATURES INCLUDE:

- Wide, full-access seat (36.75" / 93.4 cm)
- Pinch resistant flip-up L-Arms provides a safer experience for patients and caregivers.
- Universal arm design enables easy switchover to left or right-side use as your space dictates. Conversion tools included.
- Rigidly tested to an 800 lb, total load capacity in any combination of occupants – Individual, parent/child, or patient/caregiver – under normal use.
- Designed and tested to BIFMA\* standards.



*Pinch Resistant Flip-Up L-Arms*



# HARMONY BDC

## PHLEBOTOMY COLLECTION

### ANATOMY OF THE HARMONY BDC EXTRA WIDE BLOOD-DRAW CHAIR

Adjustable arm extension 17"-21" (43.1cm - 53.3cm) allows for extension of arm rest for caregiver access

Comfortable 2.5" (6.4 cm) seat & back rest

Padded L-Arms flip out of the way for easy patient access (optional)

Pinch resistant safety cover

**IMPORTANT:**  
*L-Arms are not intended to restrain a fainting or falling occupant. Never leave occupant unattended.*

**OVER 50 DESIGNER COLORS AVAILABLE**

Quick release handle adjusts arm rest height range from 4.8"-12.3" for patient and caregiver comfort

Weight distributing dual-reinforced steel center support

Footrest (optional)

Durable, powder-coated, welded tubular steel frame

6 high-impact adjustable leveling feet: stable on practically any surface



Standard Short Arm



Optional universal cabinet with table top. 6.25" x 10" deep cabinet features high impact storage drawer with spill-resistant top and rugged aluminum sides for lasting durability. Cabinet is convertible to left or right-side use as space dictates.



## SPECIFICATIONS

HARMONY BDC MODELS	2580	2580-CL	2580-CS	2580-DL	2580-SL
DESCRIPTIONS	DUAL SHORT ARMS	L-ARM - CABINET	SHORT ARM - CABINET	DUAL L-ARMS	SHORT ARM - L-ARM
Seat Height	20" (50.8 cm)	20" (50.8 cm)	20" (50.8 cm)	20" (50.8 cm)	20" (50.8 cm)
Seat Width	36.75" (93.4 cm)	36.75" (93.4 cm)	36.75" (93.4 cm)	36.75" (93.4 cm)	36.75" (93.4 cm)
Seat Depth	18.5" (47 cm)	18.5" (47 cm)	18.5" (47 cm)	18.5" (47 cm)	18.5" (47 cm)
Overall Height	35.75" (90.81 cm)	35.75" (90.81 cm)	35.75" (90.81 cm)	35.75" (90.81 cm)	35.75" (90.81 cm)
Overall Width	44" (111.8 cm)	51" (129.5 cm)	50.5" (128.3 cm)	46" (116.8 cm)	45" (114.3 cm)
Back Height above Seat	15.75" (38.1 cm)	15.75" (38.1 cm)	15.75" (38.1 cm)	15.75" (38.1 cm)	15.75" (38.1 cm)
Overall Length	24" (61 cm)	29" (73.7 cm)	24" (61 cm)	29" (73.7 cm)	29" (73.7 cm)
Adjusted Arm Range	24.5" - 32.25" (62.2 cm - 81.9 cm)	24.5" - 32.25" (62.2 cm - 81.9 cm)	24.5" - 32.25" (62.2 cm - 81.9 cm)	24.5" - 32.25" (62.2 cm - 81.9 cm)	24.5" - 32.25" (62.2 cm - 81.9 cm)
Weight	77 lbs (34.9 kg)	97 lbs (44 kg)	87 lbs (39.5 kg)	97 lbs (44 kg)	87 lbs (39.5 kg)
Weight Capacity	800 lbs (362.9 kg)	800 lbs (362.9 kg)	800 lbs (362.9 kg)	800 lbs (362.9 kg)	800 lbs (362.9 kg)

OPTIONAL FOOTREST (FR): 4.375" (11.11 cm) high x 31" (78.74 cm) wide; Adds 3.5" (8.89 cm) to overall depth and 6 lbs (2.72 kg) to overall weight  
Weight shown above is the overall weight of the chair and excludes packaging.

\* BIFMA develops, maintains, and publishes safety and performance standards for furniture products. Harmony BDC designed and tested to BIFMA X5.4-2012.

## ADDITIONAL CHAIRS IN THE PHLEBOTOMY COLLECTION



Additional options available

Spirit BDC	Solace BDC	Unity BDC
SHOWN: MODEL #2570	SHOWN: MODEL #2573	SHOWN: MODEL #2575

**WARRANTY** Contact Winco or visit our website for warranty details. **VINYL** Padded models available in over 50 vinyl colors and with TB133 Flammability option, visit our website or contact us for swatches. **NOTE** Winco reserves the right to change specifications without notice.



## CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:  
800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com





## HEALTHCARE ESSENTIALS

---

Task Seating, Treatment  
Tables & Privacy Screens



# TREATMENT TABLES

SOLUTIONS THAT FIT YOUR PRACTICE

## KEY FEATURES

- WEIGHT CAPACITY: 400 lbs (181.4 kg)
- PATIENT SURFACE (LxW): 72" (183 cm) X 28" (71 cm)
- DURABILITY: Powder-coated, steel frames
- STABILITY: Adjustable leveling feet
- FUNCTIONALITY: All models include a paper holder and are available with or without a shelf (SH)



## FLAT TOP

8500

Our most popular model provides versatility and strength. Standard table height is 30". Specialty heights (19", 25" or 36") can also be specified.



## RECOVERY COUCH

8550

Perfect for pediatric or geriatric settings. This models contoured headrest and 2.5" medium density foam provide exceptional patient comfort. Standard table height is 19", which meets ADA (Americans with Disabilities Act) maximum height guidelines.



## ADJUSTABLE BACK

8570

Features every practice needs! The gas-spring operated back features an adjustment range of 6° - 57° via a one-touch control. Standard table height is 30". Specialty heights (25" or 36") may also be specified.

ENHANCING THE PATIENT & CAREGIVER EXPERIENCE

## SPECIFICATIONS: TREATMENT TABLES

	8500 FLAT-TOP	8550 RECOVERY	8570 ADJUSTABLE BACK
HEIGHT	STANDARD: 30" (76 cm) OPTIONAL: 36" (91 cm) 25" (63 cm) 19" (48 cm)	19" (48 cm)	STANDARD: 30" (76 cm) OPTIONAL: 36" (91 cm) 25" (63 cm)
WIDTH	28" (71 cm)		
LENGTH	72" (183 cm)		
WEIGHT	STANDARD HEIGHT: 74 lb (33.6 kg) OPTIONAL HEIGHTS: [36"] 76 lb (34.5 kg) [25"] 72 lb (32.7 kg) [19"] 70 lb (31.8 kg)	73 lb (33.1 kg)	STANDARD HEIGHT: 74 lb (33.6 kg) OPTIONAL HEIGHTS: [36"] 76 lb (34.5 kg) [25"] 72 lb (32.7 kg)
SHELF	OPTIONAL *ADD 25 lb (11.34 kg) to WEIGHT		
WEIGHT CAPACITY	400 lb (181.4 kg)		



**WARRANTY** Contact Winco or visit our website for complete warranty details.

**VINYL COLORS** Visit our website to view available vinyl colors or contact our Customer Care team to review color options and request swatches.

**CERTIFICATIONS** 1) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture  
2) NR or TB options are available for upholstered items. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for upholstered furniture.



## CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:  
**800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com**



Privess™ patient privacy screens provide a broad range of flexible solutions to meet a wide spectrum of requirements for healthcare professionals, military procurement, international aid, disaster relief coordinators, and educational institutions worldwide.

- Cost-effective selections: innovative, modular options to attractive, designer models
- Permanently mounted, stationary, portable and transportable configurations available.
- Sturdy steel or lightweight aluminum construction options with smooth, easy-clean panels.

Privess™ represents our ongoing dedication to privacy: enabling clinicians to confidently focus on their patients. For every conceivable need, budget or profession, Privess™ has a custom tailored answer for your patient privacy needs.

# PRIVACY SOLUTIONS FOR A RANGE OF SPACIAL AND SITUATIONAL REQUIREMENTS



Model 3130 shown

## PRIVESS™ BASIC

### 3-PANEL PRIVACY SCREEN | 3130 / 3139

For healthcare, disaster relief and education professionals who require a sturdy metal framed screen that can easily be moved to protect patient privacy. The Privess™ Basic is a cost-effective, durable 3-panel screen that features:

- Casters for mobility or crutch-tips for stationary use
- Choice of smooth white vinyl (3130) or optional Sure-Chek® colors (3139) for easier cleaning
- Durable steel 3-panel design

Unlike similar style aluminum models the Privess™ Basic 3-Panel Privacy screen offers portable, high-quality steel framed privacy protection in a value priced package. Imported (Taiwan).



Model 3170-17 shown

## PRIVESS™ ELITE

### DESIGNER PRIVACY SCREEN | 3170

For quality and style-minded users who need a larger attractive patient privacy solution that matches their décor. Privess™ Elite is a high quality screen featuring:

- Casters for mobility or crutch-tips for stationary use
- Quality healthgrade vinyl in your choice of 50+ colors to enhance any décor
- Durable steel construction

Privess™ Elite showcases significantly upgraded durability and flexibility unlike imported lightweight models with a thin sheet of vinyl. The perfect solution for healthcare institutions, surgery centers, physician offices, and clinics where stylish substantial privacy screens are desired. Made in USA.



Model 3730 shown

## PRIVESS™ MODULAR

### 3 AND 4 PANEL CONNECTING PRIVACY SCREENS | 3730 / 3739 / 3740 / 3749

Value-minded healthcare professionals will appreciate these flexible patient privacy solutions. Privess™ Modular screens feature:

- Infinite connectivity to create the optimal size and shape for your specific needs.
- Choice of smooth white vinyl (3730 & 3740), Sure-Chek® colors (3739 & 3749) or Clear-Vue (3739-9K / 3749-9K).
- Sturdy aluminum frame and crutch tips for extra stability.

Privess™ Modular screens can swivel up to 350° when snapped-together, providing a myriad of configuration options. The Privess™ Modular's ability to fold flat coupled with its lightweight, flexible design make it ideal for rapid deployment and temporary use situations. Made in USA.





**PRIVESS™**

ESSENTIALS COLLECTION

## ANATOMY OF PRIVESS™ PRIVACY SCREENS

### PRIVESS™ BASIC AND ELITE



**Durable,  
Powder-Coated Steel Frame**


**Smooth, Easy-Clean Panels**

**Crutch Tips  
and Casters Included**

*Shown: Privess Elite™ Model 3170-17*

### PRIVESS™ MODULAR

**Sturdy,  
Aluminum Frame**



**Flat, Easy  
Clean Panels**



**Modular  
design**

**Innovative clips  
easily connect  
multiple models**

*Shown: Privess™ Modular Model 3749-94*

*Shown (clockwise):  
Privess™ Modular Models 3739-94, 3739-95 and 3739-91*

**CUSTOMIZABLE FOR DESIGN AND FUNCTION:** All models available in your choice of standard White vinyl or five optional Sure-Check® colors. Elite model also available in 50+ health-grade designer vinyl colors. Modular Models also available with Clear-Vue option.

## SPECIFICATIONS

MODEL	PRIVESS™ BASIC 3130 / 3139	PRIVESS™ ELITE 3170	PRIVESS™ MODULAR 3730 / 3739   3740 / 3749
BASE SIZE	20 in. x 27 in. (50.8 cm x 68.9 cm)	23 in. x 31 in. (58.4 cm x 78.7 cm)	15 in. x 40 in. (5.9 cm x 15.7 cm) 15 in. x 70 in. (5.9 cm x 27.6 cm)
OVERALL HEIGHT	70 in. (177.8 cm) *with casters	70 in. (177.8 cm)	69 in. (175.3 cm)
OVERALL WIDTH	79 in. (200.7 cm)	84 in. (213.4 cm)	60 in. (152.4 cm)   80 in. (203.2 cm)
SINGLE PANEL LENGTH	49 in. (124.5 cm)	55 in. (139.7 cm)	46 in. (116.8 cm)
SINGLE PANEL WIDTH	27 in. (68.6 cm)	28 in. (71.1 cm)	20 in. (50.8 cm)
WEIGHT	25 lb. (11.3 kg)	45.5 lb. (20.6 kg)	9 lb. (4.1 kg)   12 lb. (5.4 kg)

All specifications subject to change without notice, please contact our customer care team to verify any critical dimensions prior to ordering.

## FEATURES & OPTIONS

MODEL	PRIVESS™ BASIC 3130 / 3139	PRIVESS™ ELITE 3170	PRIVESS™ MODULAR 3730 / 3739   3740 / 3749
TRANSPORTABLE	Yes	Yes	Yes
MODULAR	No	No	Yes
TELESCOPING	No	No	No
FRAME	Steel	Steel	Aluminum
FEET	Crutch Tips & Casters	Crutch Tips & Casters	Crutch Tips
COLOR OPTIONS	Standard White or Sure-Chek® Colors	Standard White, Sure-Chek® Colors, or Designer Vinyl Colors	Standard White, Sure-Chek® Colors, or Clear-View
ORIGIN	Imported (Taiwan)	USA	USA

**WARRANTY** Contact Winco or visit our website for warranty details.

**COLORS** Contact Winco or visit our website for all available color selections.

**NOTE** Privess™ Privacy Screens are not intended for wet environments. Privess™ Privacy Screens are not recommended for outdoor use.

*Sure-Chek® is a registered trademark of Herculite Products Inc.*



## CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:  
**800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com**

# PRIVESS™ CLEAR-VUE

## MODULAR SCREENS



### KEY FEATURES

Privess™ Clear-Vue screens provide a see-through protective barrier solution with flexible applications.


- Cost-effective
- Lightweight, portable and configurable
- Easily cleaned and disinfected
- Compact storage requirements

Privess™ Clear-Vue provides an additional layer of protection from airborne contaminants. The Clear-Vue feature is available as an option for 3 or 4 panel Privess™ Modular Screens.

The extreme flexibility makes Privess™ Clear-Vue a perfect solution for podiums and lecterns, workstations, student desks or for queuing situations where individuals may not be able to be adequately distanced and security concerns demand visibility.

Privess™ Modular Screens feature innovative clips that allow multiple screens to easily connect to one another to create unique shape and size solutions to fit any space or need.

## SPECIFICATIONS: PRIVESS™ MODULAR SCREENS



PRIVESS™ MODULAR	3 PANEL	4 PANEL
<b>BASE SIZE</b>	15 in. x 40 in. (5.9 cm x 15.7 cm)	15 in. x 70 in. (5.9 cm x 27.6 cm)
<b>OVERALL HEIGHT</b>	69 in. (175.3 cm)	
<b>OVERALL WIDTH</b>	60 in. (152.4 cm)	80 in. (203.2 cm)
<b>SINGLE PANEL LENGTH</b>	46 in. (116.8 cm)	
<b>SINGLE PANEL WIDTH</b>	20 in. (50.8 cm)	
<b>WEIGHT</b>	9 lb. (4.1 kg)	12 lb. (5.4 kg)
<b>FRAME</b>	Aluminum	
<b>ORIGIN</b>	Made in USA	

### CONNECTIVITY

Easily connect screens to build larger screens in varying shapes and sizes.



*Shown using colored screens to better illustrate connectivity options.*

**WARRANTY** Contact a Winco representative or visit [WincoMfg.com](http://WincoMfg.com) for warranty details.

**NOTE** Privess™ Clear-Vue Screens are not intended for wet environments and should never be used as shower curtains. Privess™ Clear-Vue Screens are not recommended for outdoor use.



### CONTACT WINCO

To learn more about the Privess™ Clear-Vue option:

**800-237-3377 | [WincoMfg.com](http://WincoMfg.com) | 5516 SW 1st Lane, Ocala, FL 34474 | [Info@WincoMfg.com](mailto:Info@WincoMfg.com)**



## TASK SEATING SOLUTIONS TO MATCH YOUR NEEDS

### KEY FEATURES

- WEIGHT CAPACITY: 250 lbs (113 kg)
- SEAT DIAMETER: 16.5" (41.9 cm)
- ADJUSTABLE SEAT-HEIGHT: 17" - 22.25" (Standard & Deluxe models) 23.5" - 33.75" (Tall Models)
- PADDED UPHOLSTERY: Available in your choice of 50+ vinyl colors
- 2" TWIN WHEEL CASTERS or HIGH RESISTANCE (HR) Casters



### TASK STOOLS

**STANDARD: 4300 / 4308 | DELUXE: 4400/4408**

Multi-purpose gas-lift task stools. Standard models feature a high-density plastic base. Deluxe models feature a polished aluminum base. Upgrade any of the above models with an optional foot ring (FR).



### TASK CHAIRS

**STANDARD: 4350 / 4358 | DELUXE: 4450 / 4458 | TALL: 4380 / 4388**

Versatile caregiver seating featuring "one-touch" seat height adjustment and height-adjustable seat backs. The Standard and Tall models each include a high-density plastic base. The Deluxe models include a polished aluminum base. Our Tall model includes provides an additional 11.5" in seat height and includes a foot ring. A foot ring (FR) can also be added to any of the standard or deluxe models.



### HIGH RESISTANCE CASTERS

HR Casters allow free chair movement when weighted or occupied, but resist unnecessary chair movement when unweighted or when exiting the chair. \*HR Casters may not release and roll freely if occupant weighs less than 80lbs.

Copyright © 2021 Winco Mfg., LLC



## CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:  
**800-237-3377 | [WincoMfg.com](http://WincoMfg.com) | 5516 SW 1st Lane, Ocala, FL 34474 | [Info@WincoMfg.com](mailto:Info@WincoMfg.com)**