

## Stryker 1115 Prime Big Wheel Stretcher -Refurbished



\*Photo is stock photo and may include options not offered

The Stryker 1115 Prime Series Stretchers are extraordinarily mobile and easy to use. Big Wheel steering helps make transport safe and efficient while reducing the physical strain on the clinician.

Features include: Dual-side mounted foot controls with uni-lower pedal, dual end brake/steer controls, glideaway siderails with integrated push/pull handles, pneumatic backrest.

The Stryker 1115 Prime stretcher comes completely refurbished with a comprehensive warranty and a new pad.

- Glideaway Side Rails
- Central Brake and Steer
- $\pm$  Height,  $\pm$  Trendelenburg
- Pneumatic Backrest
- O2 Tank Storage
- 700 lbs Weight Capacity
- Optional IV Pole
- Includes New 4" Pad
- Completely Refurbished

### Specifications:

- Overall Length: 85"
- Overall Width: 34"
- Patient Surface Width: 26"
- Height Range: 20.75"– 34"
- Weight Capacity: 700 lbs
- Fowler Angle: 0° - 90°
- Trend/Rev: +16° to -17°
- Wheel Diameter: 8"
- Pad Thickness: 4"

\*Please note that head rest is tapered, but does not articulate.  
\*Optional Knee Gatch and Side Brake/Steer controls are available.