

# PlumeVac™ Smoke Evacuator

The portable PlumeVac Smoke Evacuator gives you the flexibility to put smoke plume evacuation where you need it – on or off the cart.

PlumeVac Smoke Evacuator AEVV120



- 1 Motor Speed Control**  
Three motor speeds
- 2 Filter Life Indicator**  
Provides a visual indication of filter life (max filter life 35 hours)
- 3 PlumeGuard™ Smoke Filter**  
Four stage ULPA filter captures surgical smoke and accommodates 1/4" (6 mm), 3/8" (9.5 mm) and 7/8" (22 mm) tubing diameters
- 4 Automatic Activation Device**  
Automatically activates smoke evacuator while ESU is in use



Automatic Activation Device AEAAD01

