



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Biodex 840 Imaging Table – Refurbished



Stock Photo – Options May

- Stable, accessible, Precise style, Quiet, and Vibration-Free Positioning
- Carbon Fiber tabletop
- Stainless Steel Finish
- New Mattress
- Standard OR Accessory Rails
- Warranty Included

Specifications:

Biodex 840 Imaging Table – Completely Refurbished. The Biodex 840 is engineered for image-guided fluoroscopic procedures requiring stability, precision, and unobstructed access. This fully refurbished unit supports a 500 lb patient weight capacity.

Its large radiolucent area, free of cross members, allows full fluoroscopic visualization and unrestricted C-Arm positioning. The free-float X-Y tabletop delivers quiet, vibration-free movement.

Motorized controls include height adjustment, isocentric lateral roll, and Trendelenburg motion. This refurbished Biodex 840 comes backed by a 6-month warranty.

- Tabletop Dimensions: 101" L × 24" W
- Overall Dimensions: 104" L × 26" W × 43.5" H (max)
- Height Range: 33.5" – 43.5"
- Tabletop X Motion (head-to-toe float): 35"
- Tabletop Y Motion (side-to-side float): 10"
- Trendelenburg: 0° – 20°
- Isocentric Lateral Roll: 0° – ±15°
- Patient Capacity: 500 lbs
- Unit Weight: 641 lbs
- Shipping Weight: ~865 lbs
- Mattress: New
- Warranty: 6 months
- Compliance: CDRH Registered

***Brochure and images are a general overview of the product. Actual product options may vary.*

VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReQuip, Inc.