



STRETCHER

Designed for Outstanding Mobility & Safety

A feature-rich standard configuration to meet your needs and your budget.



Fold-away push handles
improve maneuverability



Sliding Backrest
reduces patient migration resulting in less boosting



Head of bed angle indicator
aids with airway clearance head elevation compliance.



Knee gatch
provides the patient added comfort. Easy to use manual crank.



Fifth wheel
allows optimal handling and control during transport.

**Three-year
Limited Product
Warranty**

**700 lbs Safe
Working Load**

STRETCHER



Examination and Recovery Configuration

Includes Knee Gatch
Recovery Position

Transport Configuration

Stationary foot-end
patient platform

Standard Features

- 700 lbs weigh capacity
- 28-inch width patient deck
- Fifth wheel
- Knee Gatch /Recovery Position (on Examination and Recover configuration only *)
- Fold-away push handles on head-end
- Sliding backrest
- Permanently mounted IV pole at head-end, 2 hooks
- Angle indicators for head of bed and trend
- 8-inch integral casters with wheel covers
- Pneumatic assisted backrest
- Trendelenburg and reverse trendelenburg
- Central-locking brake system
- Brake and steer pedals on head and foot ends
- Hydraulic elevation control pedals on both sides
- Fold down siderails
- 4 integrated IV sockets
- Drainage bag holders
- Bumper wheels on each corner
- Patient storage tray on base
- Oxygen tank holder on base
- 4" Comfort Series Pad

Specifications

- Weight Capacity 700 lb
- Patient surface width 28"
- Patient surface length 77"
- Knee gatch 0° - 45
- Backrest 0° - 90°
- Height range 23.25" - 34"
- Trendelenburg 17°
- Reverse trendelenburg -17°
- Caster diameter 8"
- Overall length 84"
- Overall width 36"

Optional Features

- Deck Mounted Oxygen Tank Holder
- Fold-away push handles on foot-end
- Side Rail Pads
- Chart holder
- Removable IV Pole, 2 stage and 2 hooks

Warranty

Three-year limited product warranty

MCM9010	Premium Transport/Recovery Stretcher, Knee Gatch, IV Pole, Base Oxygen Tank Holder, 700 lbs. weight capacity
MCM9022	Foldable Push Handles on Foot End
MCM9023	Removable IV Pole, adjustable height, 2 Hooks
MCM9024	Chart Holder
MCM9025	Side Rail Pads
MCM9026	Deck Mounted Oxygen Tank Holder