

MI 1000 LED

Bovie®

MI 1000 LED

THE NEXT GENERATION
IN SURGERY LIGHTING



COMPELLING

The MI 1000 light is part of the "cutting-edge" MI LED SERIES by Medical Illumination International

COOL

The MI 1000 light emits very low heat

BRIGHT

The MI 1000 light is output rated at 100,000 lux (9,300 foot candles) at 4,300 degrees Kelvin color temperature

EFFICIENT

"GREEN TECHNOLOGY"...The MI 1000 light uses significantly less electric consumption than standard halogen lights

RELIABILITY

Lighting equipment (excluding replacement parts) has a three (3) yr. limited warranty. Long-life LED's, rated at 50,000 hrs.

VERSATILE

Available in a variety of mounting options

INNOVATIVE


Five-stage dimming and on/off switch controlled from sterilizable handle

CONVENIENT

No replacement bulbs or fuses necessary...lower overall maintenance costs

SURGICAL LIGHTING

MI 1000 LED
INCLUDES



OTHER MI 1000 LED MODELS
AVAILABLE



For more product information, contact your Bovie sales representative.

MI 1000 LED TECHNICAL
SPECIFICATIONS

Dimensions:		Rotations:	
Head	18" (457 mm) Dia. X 4" High (102 mm)	Articulating Arm Vertical Movement	+/- 40 Degrees travel
Arm	57.0" (1448 mm) Long	Articulating Arm Horizontal Movement	Approximately 540 Degrees
Arm (Floor)	38.5" (978 mm) Long	Articulating Arm/Yoke Interface	Approximately 540 Degrees
Floor Stand	74.5" (1892 mm) High	Yoke/Lamp Head Interface	Approximately 190 Degrees
Weight:		Data:	
Head	Approximately 13.5 lbs. (6.1 kg)	Light Source	LED
Arm	Approximately 12.5 lbs. (5.7 kg)	LED Life in Hours	50,000 (Average)
Arm (Floor)	Approximately 11.5 lbs. (5.2 kg)	Kelvin Color Temperature	4300
Wall Bracket Assy.	Approximately 3.0 lbs. (1.4 kg)	Lux	100,000
Floor Stand Assy.	Approximately 43.0 lbs. (19.5 kg)	Power Consumption	60 Watts
Single Ceiling Mount	Approximately 14.5 lbs. (6.6 kg)	Input Voltage	100 - 240 VAC 50/60 Hz
Dual Ceiling Mount	Approximately 16.0 lbs. (7.3 kg)		

Made in the USA. The MI LED Series is currently only available for sale in the USA.