



stryker®

Image Quality.
Ultimate Access.
Versatility.

Stryker Prime X™ X-ray Stretcher

Prime Series® Stretchers

Stryker stretchers have been the industry benchmark for decades, and Stryker's Prime Series Stretchers continues the tradition.

The Stryker Prime Solution, Now with X-Ray Capability

Advanced mobility, safety and efficiency combine with the Clearview Technology platform.

The Prime X X-Ray Stretcher is the next generation of Prime Series Stretchers, loaded with all of the features you already love — now, with X-ray capability. Stryker's exclusive Clearview Technology platform provides similar image quality to that of a stationary X-ray table, on the Prime Series Stretcher mobile platform.

- 360 degree access for cassette placement
- 52% larger imaging window
- 700 lb. weight capacity
- Integrated scale
- Power-washable design

The results? Increased safety; incredible efficiency.

Electric Lift

With electric lift, simply touch the pedal to raise the litter. Available on Fifth Wheel and Big Wheel Advanced Mobility models.





The Prime X X-ray Stretcher offers the safety and efficiency of the Prime Series Stretcher, now with the Clearview Technology platform.






that provides controlled
mobility which

NOTE: Options not available in all countries. Please contact your local
country representative to learn of availability in your area.

The Prime X X-ray Stretcher allows you to **adjust the cassette to the patient** rather than the patient to the cassette.

The unique dual-deck design allows positioning of cassettes at any point under the patient, from head to foot and side-to-side. Open architecture accommodates virtually any size of cassette.



Easy to Raise Foot Section offers enhanced cleanability through convenient access to the cassette platform.



Open Cassette Platform and a full 360° of access allows positioning of cassettes from head end to foot end.



Beveled Edges and Visual Alignment Guides make positioning the cassette to the patient simple.



Mobile X-ray platform helps reduce patient transfers, offering a versatile solution throughout the hospital.



Reconstructive

Hips
Knees
Trauma & Extremities
Joint Preservation
Orthobiologics

Medical & Surgical

Power Tools & Surgical Accessories
Image Guided Navigation
Endoscopy & Arthroscopy
Integrated Communications
Beds, Stretchers & EMS
Sustainability Solutions

Neurotechnology & Spine

Craniomaxillofacial
Interventional Spine
Neurosurgical, Spine & ENT
Neurovascular
Spinal Implants

3800 E. Centre Ave.
Portage, MI 49002 U.S.A.
t: 269 329 2100 f: 269 329 2311
toll free: 800 327 0770

www.stryker.com

The information presented in this brochure is intended to demonstrate Stryker products. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker Account Manager if you have questions about the availability of Stryker products in your area.

Products referenced with [™] designation are trademarks of Stryker.
Products referenced with [®] designation are registered trademarks of Stryker.

Mkt Lit-625 13 DEC 2011 Rev B.0
MHD/PS **3M** 1714, SPH255

Copyright © 2014 Stryker
Printed in U.S.A.