



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Biodex 240 MRI Stretcher - Refurbished



Stock Photo – Options May Vary

- Heavy Duty MRI Stretcher
- All Aluminum, Non-magnetic
- Includes New Mattress
- Retractable Side Rails
- Perfect for Patient Transport
- Includes Patient Restraint Straps
- Completely Refurbished

The Biodex 240 MRI Stretcher is a heavy-duty, all aluminum non-magnetic stretcher and is an ideal unit for MRI patient transport.

Although light in weight for easy, controlled maneuverability, the MRI stretcher accommodates up to a 500 lb patient safely and comfortably.

Standard features include: retractable guard rails, patient restraint straps, and easy to roll 5" locking, swivel casters.

The Biodex 240 MRI stretcher comes completely refurbished with comprehensive warranty and a new pad.

Specifications:

- Tabletop: 72" x 24"
- Overall: 74.5" x 27.75" x 34"
- Wheels: 5" locking casters
- Finish: Powder coat
- Shipping Weight: 200 lbs
- Patient Capacity: 500 lb
- Weight: 85lb

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReQuip, Inc.