



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Hausted 4160 Pediatric Crib - Refurbished



Stock Photo – Options May Vary

The Hausted 4160 Pediatric Stretcher/Crib is ideal for patients from infant to 6 years old. Three position vertical side rails lock in place where needed and the dual pedestal design provides maximum stability for the patient.

The 4160 is height adjustable and also can be placed into either trendelenburg or reverse trendelenburg. The backrest is manually operated.

The Hausted 4160 Pediatric stretcher/crib comes completely refurbished with a comprehensive warranty and a new pad.

- Includes New Pad
- Hydraulic Height Adjustment
- Manual Backrest
- Plastic Base with Storage
- Includes O2 Tank Well
- Full Length Side Rails
- Central Brake and Steer
- Weight Capacity 325 lbs
- Completely Refurbished

Specifications:

- Minimum Height: 32.5"
- Maximum Height: 39"
- Overall Length: 59.25"
- Overall Width: 29.75"
- Patient Surface: 25.25" x 56"
- Side Rails: 17.75" x 56"
- End Rails: 17.75" x 23"
- Caster Diameter: 8"
- Trendelenburg : 10°
- Fowler Angles: 0-45°
- Weight Capacity: 325 lbs

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReQuip, Inc.