



www.venturemedical.com

6008 Bonacker Drive

Tampa, FL 33610

800-627-3215

813-740-8640

Fax: 800-466-1251

sales@venturemedical.com

Hausted Horizon 462 Airglide Stretcher - Refurbished



Stock Photo – Options May Vary

- Airglide Tuck Under Side Rails
- Brake and Steer Functions
- $\pm$ Height, $\pm$ Trendelenburg
- Gas Spring Operated Fowler
- 500 lbs Weight Capacity
- Dual Side Foot pedal Controls
- Includes New 4" Pad
- Completely Refurbished

The Hausted Horizon Series AirGlide® Stretcher features include Chrome Plated, Full Length, Tuck Under AirGlide® Patient Safety Rails with Dual Pneumatic Assists, Dual Side Foot Pedal Hydraulic Controls for Height and Trendelenburg Adjustment, and Pneumatic Assist Backrest. Brake and Steer Functions can be actuated From Any Corner. The ABS Plastic Base Cover has a storage area and O2 Tank Holder.

The 500 lbs weight capacity can handle large patients with ease.

The Hausted Horizon 462 comes completely refurbished with a comprehensive warranty and a new pad.

Specifications:

- Overall Length: 84"
- Overall Width: 31.5"
- Overall Width (Wide Model): 34.75"
- Patient Surface: 28"W x 76"L
- Patient Surface (Wide Model): 32"W x 77"L
- Height Range: 22.25" – 34.75"
- Weight Capacity: 500 lbs
- Fowler Angle: 0° - 90°
- Trend/Rev Trend: $\pm$ 16°
- Side Rails: 15" x 77"
- Wheel Diameter: 8"
- Pad Thickness: 4"

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReQuip, Inc.