



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Hausted 493 Fluoro-Trak Stretcher - Refurbished



Stock Photo – Options May Vary

Hausted 493 Fluoro-Trak Stretcher: Refurbished. The Hausted Fluoro-Trak Fluoroscopy Capable Stretcher is ideal for fluoroscopy, endoscopy, and other procedures. The Fluoro-Trak has a C-Arm clearance of 24" and has a patient weight capacity of 500 lbs. The 493 features Sliding Top that travels a total of 12" for positioning during procedures and radiolucent patient surface.

Additional 493 Fluoro-Trak Features: The 493 features Adjustable Height, Trendelenburg, and Pneumatic Assist Backrest. Travels easily on 8" Casters with Four Brake/Steer Levers located at each corner for easy implementation and steers from all four corners. Hydraulic controls: foot pedals on both sides and foot end IV and accessory wells, one at each corner. Stainless Steel 6-leg, 3 position retractor rails with built in Push/Steer Handles, and top caps to eliminate patient transfer gaps. IV and accessory wells: One IV and one Accessory Well at each corner

- C-Arm capable
- Radiolucent Patient Surface
- Pneumatic-assist Backrest
- Adjustable Height
- Trendelenburg/Reverse Trendelenburg
- Sliding Top
- Hydraulic controls: foot pedals on both sides
- 8" Casters
- Foot End IV and Accessory Wells
- 3" Mattress Pad - New
- Weight Capacity: 500 lbs.

Specifications:

- Fluoroscopy Capable Stretchers: 0-80°/ 29.75", 15", 36"
- Patient Surface: Radiolucent
- C-Arm clearance: 24"
- Sliding top: 12" toward foot end and 6" toward head end
- Trend/Reverse Trend: 15°
- Pneumatic-assist backrest: 0-80°
- Height: 23" - 36"
- Overall Dimensions: 29.75"W X 80.25"L
- Hydraulic controls: foot pedals on both sides and foot end
- Brake/Steer Levers: 4 - one at each corner
- Casters: 8"
- Side rails Dimensions: 15" x 64"

***Brochure and images are a general overview of the product. Actual product options may vary.*

VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReQuip, Inc.