



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Hill-Rom 8050 OB-GYN Stretcher - Refurbished



Stock Photo – Options May Vary

- Gas Spring Assist Backrest
- Collapsible Side Rails
- Central Brake and Steer
- $\pm$ Height, $\pm$ Trendelenburg
- Built-in O2 Tank Holder
- 500 lbs Weight Capacity
- Includes New 4" Pad
- Completely Refurbished

The Hill-Rom 8050 OB-GYN stretcher offers outstanding patient positioning, quickly and easily. The stretcher comes with integrated foot supports and a sliding top. The foot section retracts out of the way during procedures, leaving the footrests exposed.

Standard features include lift assist backrest (Fowler), Trendelenburg, reverse Trendelenburg, central brake and steer, height adjustment, and 500 lbs weight capacity.

The Hill-Rom 8050 OB-GYN stretcher comes completely refurbished with a comprehensive warranty and a new pad.

Specifications:

- Overall Length: 80.5"
- Overall Width: 35.5"
- Patient Surface Width: 30"
- Height Range: 24 $\frac{1}{4}$ " - 38"
- Weight Capacity: 500 lbs
- Pad Thickness: 4"
- Side Rail Length: 37"
- Caster Size: 8"

***Brochure and images are a general overview of the product. Actual product options may vary.*

VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReQuip, Inc.