



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Luvis E100 LED Exam Light - New



Stock Photo – Options May Vary

- Excellent Illumination
- Ergonomic & Friendly Design
- No UV Emission & Cool Light
- Excellent Depth of Illumination
- 3 Selectable Color Temperature Options
- Intuitive and Easy Way of Control
- Intensity: 60,000 lux
- Color Rendering Index: 95
- Lightweight Suspension System
- Various Installation Options
 - Ceiling, Mobile, Wall

The Luvis E100 LED Exam Light is ideal for medical examinations in smaller spaces. The E100 delivers value through longer life, lower energy consumption and will save the cost of hassled maintenance.

The E100 LED Exam Light features three selectable color temperature options (3,500K, 4,000K, and 4,500K).

The Luvis E100 Exam Light can come in ceiling, mobile and wall configurations and comes with a 2-year mfg. warranty.

Specifications:

- LED's: 8
- CRI: 95
- Color Temperature:
 - 3,500K /4,000K/4,500K
- Head Size: 9.84" x 9.68"
- Focus Size: 7" (Circle)/11.8" (Circle)
- Focal Distance: 1M
- Intensity: Max. 60,000lux
- Input Power: 90-264Vac, 50/60Hz
- Power Consumption: 20W
- LED Life Span: 50,000hr

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM