



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Luvis L400 LED Dual Surgical Light - New



Stock Photo – Options May Vary

- Excellent Illumination
- Open Center Cavity
- Ergonomic & Modern Design
- Detachable & Autoclavable Main Handle
- No UV Emission & Cool Light
- Excellent Depth of Illumination
- 4 Step Color Temperature (selectable)
- Intuitive and Easily Controlled by User
- Rotating Handle for Light Field Size
- Intensity: 160,000 lux
- Color Rendering Index: 97
- Focusable Light Head
- Lightweight Suspension System
- Various Installation Options
 - Ceiling, Mobile, Dual Ceiling Mount

The Luvis L400 LED Surgical Light is not only an ultra-modern, full featured major surgical light, but also has an open center cavity, which promotes air flow over the surgical site. The L400 boasts exceptional quality with energy efficient operation, and the modern features you would expect for a major surgical light in today's OR environment.

The 96 LED'S in each light head produce up to 160,000 lux in a pattern that can be adjusted from 7.8" to 11.8". The Luvis Surgical Light also features a light head mounted control panel that can be synchronized together (dual config) for easy adjustment during surgery.

The Luvis L400 surgical light can come in a single, dual or mobile configuration and includes a 2-year mfg. warranty.

Specifications:

- LEDs: 96
- CRI: 97
- Color Temperature: 3.6k, 4.2k, 4.8k, 5.4k
- Head Size: 28.7"
- Focus Size: 7.8" (Circle)/11.8" (Circle)
- Focal Distance: 1M
- Intensity: Max. 160,000lux
- Input Power: 100-264Vac, 60Hz
- Power Consumption: 125W
- LED Life Span: 50,000hr

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM