



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Luvis M200 LED Surgical Light - New



Stock Photo – Options May Vary

- Excellent Illumination
- Ergonomic & Friendly Design
- Detachable & Autoclavable Main Handle
- No UV Emission & Cool Light
- Excellent Depth of Illumination
- 4,300k Color Temperature
- Intuitive and Easy Way of Control
- Rotating Central Knob
- Intensity: 120,000 lux
- Color Rendering Index: 95
- Focusable Light Head
- Lightweight Suspension System
- Various Installation Options
 - Ceiling, Mobile, Dual Ceiling Mount

The Luvis M200 LED Surgical Light is the perfect light for any procedure. The M200 features exceptional quality with energy efficient operation.

The 54 LED'S in each light head produce up to 120,000 lux in a pattern that can be adjusted from 7.8" to 11.8". The Luvis Surgical Light also features a control panel that is synchronized with sub-handle for the user to approach the control panel easily during surgery.

The Luvis M200 surgical light can come in a single, dual or mobile configuration and includes a 2 year mfg. warranty.

Specifications:

- LED's: 54
- CRI: 95
- Color Temperature: 3,800k/4,300k/4,800k
- Head Size: 19.6"
- Focus Size: 7.8" (Circle)/11.8" (Circle)
- Focal Distance: 1M
- Intensity: Max. 120,000lux
- Input Power: 100-240Vac, 60Hz
- Power Consumption: 85W
- LED Life Span: 50,000hr

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM