



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Luvis M300 Dual Ceiling Mounted LED Surgical Light - New



Stock Photo – Options May Vary

The Luvis M300 LED Surgical Light is the perfect light for any procedure. The M300 features exceptional quality with energy efficient operation.

The 84 LED'S in each light head produce up to 160,000 lux in a pattern that can be adjusted from 7.8" to 11.8". The Luvis M300 Surgical Light also features a control panel that is synchronized a sub-handle for the user to utilize the control panel easily during surgery.

The Luvis M300 surgical light can come in a single, dual or mobile configuration and includes a 2 year mfg. warranty.

- Excellent Illumination
- Excellent Dilution Effect
- Ergonomic & Friendly Design
- Detachable & Autoclavable Handle
- No UV Emission & Cool Light
- Excellent Depth of Illumination
- 4,300k Color Temperature
- Intuitive and Easy Way of Control
- Rotating Central Knob
- Intensity: 160,000 lux
- Color Rendering Index: 95
- Focusable Light Head
- Lightweight Suspension System
- Various Installation Options
 - Ceiling, Mobile, Dual Ceiling Mount

Specifications:

- LEDs: 84
- CRI: 95
- Color Temperature: 4,300k
- Head Size: 23.6"
- Focus Size: 7.8" (Circle)/11.8" (Circle)
- Focal Distance: 1M
- Intensity: Max. 160,000lux
- Input Power: 100-240Vac, 60Hz
- Power Consumption: 125W
- LED Life Span: 60,000hr

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM