



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Luvis S200 Minor LED Surgical Light Dual Ceiling Mounted - New



Stock Photo – Options May Vary

- Excellent Illumination
- Ergonomic & Friendly Design
- Detachable & Autoclavable Main Handle
- No UV Emission & Cool Light
- Excellent Depth of Illumination
- 3 Selectable Color temperature Options
- Intuitive and Easy Way of Control
- Excellent Scialytic Effect
- Rotating Central Knob
- Intensity: 100,000 lux
- Color Rendering Index: 90
- Focusable Light Head
- Lightweight Suspension System
- Various Installation Options
 - Chair, Ceiling, Mobile, Wall, Dual Ceiling Mount

The Luvis S200 LED Surgical Lights are the perfect light for minor procedures and examinations in smaller spaces. The S200 provides a shadowless effect that offers less shadow even under a 50% screened condition as well as a fan-less structure that generates no heat or noise.

The S200 LED Operating Lights feature an adjustable size pattern of 18cm pattern for posterior and a 30cm for anterior regions.

The Luvis S200 surgical light can come in a single or dual head configuration.

Specifications:

- LED's: 18
- CRI: 90
- Color Temperature:
 - 3,800K /4,300K/4,800K
- Head Size: 13.78" x 13.39"
- Focus Size: 7" (Circle)/11.8" (Circle)
- Focal Distance: 1M
- Intensity: Max. 100,000lux
- Input Power: 90-264Vac, 50/60Hz
- Power Consumption: 25W per head
- LED Life Span: 50,000hr

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReQuip, Inc.