



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Maquet Alphastar 1132.01 Surgery Table - Refurbished



Stock Photo – Options May Vary

- C-Arm Compatible Table Top
- Wired Hand Control
- 500 lbs Weight Capacity
- Automatic Electric Floor Locks
- Return-to-Level Function
- Battery Operated
- Accessory Side Rails
- Includes New Pad Set
- Completely Refurbished

The Maquet Alphastar 1132 is an ultra-modern operating table and a great choice for almost any surgical setting. The 1132 features included: wired hand control, battery operation, automatic electric floor locks, and return-to-level function

Additional functions include: $\pm$ back, $\pm$ leg, $\pm$ height, $\pm$ tilt, $\pm$ flex, trend/rev trend, and comes with North American side rails which fit most accessories.

The Maquet Alphastar 1132.01 comes completely refurbished with a comprehensive warranty and a new pad set.

Specifications:

- Overall Length: 82½"
- Overall Width: 22"
- Patient Surface Width: 20"
- Height Range: 27" - 44"
- Weight Capacity: 500 lbs
- Reverse Orientation: 300 lbs
- Trend/Rev Trend: $\pm$ 30°
- Back Rest: + 68° / - 25°
- Leg Rest: + 80° / - 105°
- Lateral Tilt: 18°
- Flex and Reflex Capable
- Power: 110-240 VAC
- Pad Thickness: 2"

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM