



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Luvis S300 LED Dual Surgical Light - New



- Excellent Illumination
- Synchronized Control Capable
- Offset Handle
- Ergonomic & Modern Design
- Detachable & Autoclavable Main Handle
- Cool Light
- Excellent Depth of Illumination
- 3 Step Color Temperature (selectable)
- Intuitive and Easily Controlled by User
- Light Field Size Adjustment (3 step)
- Intensity: 130,000 lux
- Color Rendering Index: 96
- Focusable Light Head
- Lightweight Suspension System
- Various Installation Options
 - Ceiling, Mobile, Dual Ceiling Mount

The Luvis S300 LED Surgical Light is an ultra-modern, full featured surgical light. The S300 boasts energy efficient operation, high lux output (130,000) and the modern features you would expect from a surgical light in today's OR environment.

Specifications:

The pattern can be adjusted from 5.9" to 9.85" using a touch pad on the light head controller. The Luvis Surgical Light also features a light head mounted control panel that can be synchronized together (dual config) for easy adjustment during surgery using an optional Android tablet.

The Luvis S300 surgical light can come in a single, dual or mobile configuration and includes a 2 year mfg. warranty.

- LEDs: 42
- CRI: 96
- Color Temperature: 3.8k, 4.3k, 4.8k
- Head Size: 17.2"
- Focus Size: 5.9" – 9.85" (Circle)
- Focal Distance: 1M
- Intensity: Max. 130,000lux
- Input Power: 100-240Vac, 60Hz
- LED Life Span: 60,000hr

VENTUREMEDICAL.COM

© 2007-2020 Venture Medical ReQuip, Inc.