



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Skytron 3600B UltraSlide Operating Table - Refurbished



Stock Photo – Options May Vary

The Skytron 3600B UltraSlide surgical table is a multi-purpose operating table that features the ability to perform full body imaging. This table is a great choice for vascular, neuro, urology and orthopedics.

The 3600B UltraSlide features include: a wired hand control, 23" (total) top slide capability, electro-hydraulic brake, C-Arm compatibility, and return-to-level, along with battery operation.

Additional functions include: $\pm$ back, $\pm$ leg, $\pm$ height, $\pm$ tilt, trend/rev trend, and comes with North American side rails which fit most accessories.

The Skytron 3600B UltraSlide comes completely refurbished with a comprehensive warranty and a new pad set.

- C-Arm Compatible Tabletop
- Wired Hand Control
- Battery Operation
- Return-to-Level Function
- Full Slide-Top
- 1000 lbs Weight Capacity (lift)
- Electro-Hydraulic Floor Locks
- North American Side Rails
- Includes New 3 Piece Pad

Specifications:

- Overall Length: 82"
- Overall Width: 22"
- Patient Surface Width: 19.50"
- Height Range: 24" – 42.5"
- Weight Capacity: 1000 lbs lift, 800 lbs articulated
- Trend/Rev Trend: $+ 30^\circ$ / -30°
- Back Rest: $+ 90^\circ$ / $- 40^\circ$
- Leg Rest: $+ 15^\circ$ / -100°
- Lateral Tilt: $\pm 30^\circ$
- Slide to Head: 9.5"
- Slide to Foot: 13.5"
- Power: 120 V
- Pad Thickness: 2"

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM