



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Stryker Cub FL19H Pediatric Crib Refurbished



Stock Photo – Options May Vary

The Stryker Cub (FL19H) Pediatric Crib provides a safe, comforting environment for the patient, while giving outstanding access for the clinician. The Stryker Cub also has an ultra-modern design with a warm, friendly appearance.

The crib features a retractable backrest, retractable fifth wheel for easy steering, height adjustment, trend/rev trend, O2 tank holder, and a new 4" patient pad.

All four side rails can be lowered to provide unobstructed access to the patient – 360° Access. All activation points to lower the side rails are clearly marked red and require two hand movements to operate for safety. The side rails have a counterbalance system integrated so they are easy to actuate on move. Lastly, swing-open access doors on each corner allow for quick and convenient patient care.

The Stryker Cub Pediatric Crib (FL19H) comes completely refurbished with a comprehensive warranty and a new pad.

- All-In-One Efficiency
- Fully Lowerable Side Rails
- Access Doors at Each Corner
- Counterbalanced Side Rails
- 5th Wheel Steering
- Oxygen Tank Holder
- Central Brake and Steer
- ± Height (foot operated)
- ± Backrest
- ± Trendelenburg (foot operated)
- 400 lbs Weight Capacity
- Includes New 4" Pad
- Completely Refurbished

Specifications:

- Model: FL19H
- Overall Length: 65.5"
- Overall Width: 37"
- Patient Surface Width: 30" x 58"
- Weight capacity: 400 lbs
- Height range: 33" - 40"
- Backrest: 0° to 50°
- Trend/Rev Trend: +/- 9°
- Foot Section Angle: 0°/6°/12°
- Caster Diameter: 6"

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReQuip, Inc.