



www.venturemedical.com

Venture Medical ReQuip, Inc.  
6008 Bonacker Drive  
Tampa, FL 33610  
800-627-3215 | 813-740-8640  
Fax: 800-466-1251  
sales@venturemedical.com

## Midmark/Pedigo 500 Crib Stretcher - Refurbished



*Stock Photo – Options May Vary*

The Midmark/Pedigo 500 Pediatric Stretcher/Crib, with its unique oval design, rolls easily on 5" oversize locking casters.

Standard features include manual fowler, one-handed operation, lowers and raises side rails with side rail intermediate stop position, three-position locking fowler, Trendelenburg, reverse Trendelenburg, and Raising Backrest.

The Midmark/Pedigo 500 Pediatric stretcher comes completely refurbished with a comprehensive warranty and a new pad.

- Includes New Pad
- One-Handed Side Rail Operation: Raise and Lower
- Side rail intermediate stop position
- Raising Backrest
- Trendelenburg / Fowler
- Completely Refurbished

### Specifications:

- Top Rail Height: 59 ½"
- Height with Top: 74"
- Overall Length: 60"
- Overall Width: 30 1/4"
- Spring Height: 32 1/4"
- Caster Diameter: 5"
- Side Rail Intermediate Stop Position: 13 3/4"
- Trendelenburg Angles: - 5,0,5,10 degrees
- Fowler Angles: 0,15,30,45 degrees
- Weight Capacity: 70.0 lbs

*\*\*Brochure and images are a general overview of the product.  
Actual product options may vary.*

**VENTUREMEDICAL.COM**