



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Skytron 3602 UltraSlide Surgical Table - Refurbished



Stock Photo – Options May Vary

The Skytron 3602 UltraSlide is a versatile refurbished surgical table designed for advanced imaging and positioning. It features a 23" tabletop slide that allows easy imaging access without reducing performance, maintaining a 1,000 lb lift capacity and 800 lb articulation capacity.

The table offers removable leg and back sections for imaging flexibility, with 9.5" slide-to-head and 13.5" slide-to-foot positioning to support full-body imaging and reduce the need for additional accessories.

Additional features include accessory side rails, a removable headrest, a powered kidney bridge, and an internal rechargeable battery backup.

The Skytron 3602 UltraSlide comes completely refurbished with new pads and is ideal for vascular, neuro, urology, and orthopaedic procedures.

- Wired Hand Control
- Internal Battery Backup (Rechargeable)
- Anatomically Correct Return to Level
- Anatomically Correct Beach Chair Position
- Slide Top Movement Capability
- Full Movement Range (Back, Height, +/- Trend, Leg, Tilt)
- Electrically Operated Kidney Bridge
- Self-Leveling Floor Locks
- Excellent C-Arm Access
- 1000 lbs Weight Capacity

Specifications:

- Weight Capacity: 1,000 lbs lift / 800 lbs articulated
- Tabletop Slide: 23" total (9.5" to head, 13.5" to foot)
- Height Range: 24" – 42.5"
- Imaging & Positioning: $\pm 30^\circ$ lateral tilt, 30° Trendelenburg / reverse Trendelenburg
- Back Section: 90° up / 40° down
- Leg Section: 15° up / 100° down
- Flex / Reflex: 40° / 90°
- Powered Kidney Bridge: 5.75"
- Table Dimensions: 82" L (with headrest) / 22" W overall

***Brochure and images are a general overview of the product. Actual product options may vary.*

VENTUREMEDICAL.COM