



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Skytron 6001 Surgical Table - Refurbished



Stock Photo – Options May Vary

- C-Arm Compatible Table Top
- Wired Hand Control
- 500 lbs Weight Capacity
- 180° Table Top Rotation
- Floor Locks
- Stainless Steel Base and Column
- North American Side Rails
- Swivel Casters
- Includes New Pad Set
- Completely Refurbished

The Skytron 6001 is a modern multi-purpose surgical table and a great choice for almost any surgical setting and is C-arm compatible. The 6001 features a wired hand control, hydro-electric brakes, and a 500 lbs weight capacity.

The 6001 operating table has all standard functions, ± back, ± leg, ± height, ± tilt, ± flex, trend/reverse trend and comes with North American side rails which fit most accessories.

The Skytron 6001 comes completely refurbished with a comprehensive warranty and a new pad set.

Specifications:

- Overall Length: 75"
- Overall Width: 21 ¾"
- Patient Surface Width: 20"
- Height Range: 28" – 45"
- Weight Capacity: 500 lbs
- Trend/Rev Trend: ± 30°
- Back Rest: + 90° / - 40°
- Leg Rest: + 80° / -105°
- Lateral Tilt: 30°
- Flex and Reflex Capable
- Power: 120 Volts
- Pad Thickness: 2"

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM