



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Skytron 6500 Elite Operating Table - Refurbished



Stock Photo – Options May Vary

The Skytron 6500 Elite operating table is a modern multi-purpose surgical table and a great choice for almost any surgical setting. The 6500HD features include wired hand control, 180° tabletop rotation, hydro-electric brakes, C-Arm compatibility, and a 500 lbs. weight capacity in normal position.

The Elite operating table has all standard functions ± back, ± leg, ± height, ± tilt, ± flex, trend/rev trend, and comes with North American side rails which fit most accessories.

The Skytron 6500 Elite comes completely refurbished with a comprehensive warranty and a new pad set.

- C-Arm Compatible Table Top
- Wired Hand Control
- 500 lbs Weight Capacity
- 180° Table Top Rotation
- Electro-Hydraulic Floor Locks
- Stainless Steel Base and Column
- North American Side Rails
- Swivel Casters
- Includes New 3 Piece Pad

Specifications:

- Overall Length: 76"
- Overall Width: 21 ¾"
- Patient Surface Width: 20"
- Height Range: 26" - 47"
- Weight Capacity: 500 lbs
- Trend/Rev Trend: ± 30°
- Back Rest: + 90° / - 40°
- Leg Rest: + 20° / -105°
- Lateral Tilt: 30°
- Flex and Reflex Capable
- Power: 120 Volts
- Pad Thickness: 2"

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReQuip, Inc.