



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Skytron 6700B Hercules Bariatric Surgical Table - Refurbished



210° Top Rotation

Stock Photo – Options May Vary


The Skytron 6700B Hercules operating table is a modern multi-purpose surgical table and a great choice for almost any surgical setting. The Skytron 6700B features include wired hand control, 210° tabletop rotation, hydro-electric brakes, C-Arm compatibility, and 1200 lbs. weight capacity (lift) in normal position. Even when articulated, the weight capacity is still an impressive 1000 lbs.

The 6700B also has a low height of 23", making transfers safer and easier. The leg and back section are removeable, which makes anesthesia access easier for shorter patients in lithotomy position, and more legroom for seated surgeons. Removeable leg section also creates more free space for cysto and urology procedures.

The Skytron 6700B has all standard functions $\pm$ back, $\pm$ leg, $\pm$ height, $\pm$ tilt, $\pm$ flex, trend/rev trend, and return-to-level. The 6700B is also equipped with battery operation.

This Skytron 6700B Hercules surgical table comes completely refurbished with a comprehensive warranty and new pad set.

- C-Arm Compatible Tabletop
- Battery Backup
- Wired Hand Control
- 1200 lbs Weight Capacity (lift)
- 1000 lbs Capacity (articulated)
- 210° Tabletop Rotation
- 23" Low Table Height
- Removeable Back and Leg Sections
- Auxiliary Controls on Base
- Electro-Hydraulic Floor Locks
- Stainless Steel Base and Column
- North American Side Rails
- Swivel Casters
- Includes New 3 Piece Pad



Specifications:

- Overall Length: 79"
- Overall Width: 22"
- Patient Surface Width: 19.75"
- Height Range: 23" - 44"
- Weight Capacity: 1200 lbs lift
- Trend/Rev Trend: $\pm$ 30°
- Back Rest: + 90° / - 40°
- Leg Rest: + 20° / - 100°
- Lateral Tilt: 30°
- Flex and Reflex Capable
- Power: 120 Volts
- Pad Thickness: 2"

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM