



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Skytron 6701 Hercules Surgical Table - Refurbished



Stock Photo – Options May Vary

- Full Body Imaging
- Wired Hand Control
- 1,200 lb. Weight Capacity
- 1,000 lb. Weight Capacity when articulated.
- 210° Table Top Rotation
- 23" Low Height
- Auto Beach Chair position
- North American Side Rails Fit Most Accessories
- New Pads
- Completely Refurbished
- (Lift), 800 Lbs (articulated)

The refurbished Skytron 6701 Hercules Surgical Table delivers strength and versatility for a broad range of procedures. With a 1,200 lb. weight capacity, it supports bariatric and general surgery patients while maintaining stability and precise control.

A 23" low height adjustment improves access and ergonomics, and 210° tabletop rotation allows full-body imaging without repositioning. The Self-Leveling Brake Safety System ensures secure locking during use, and North American side rails provide wide accessory compatibility.

Ideal for plastic surgery, orthopedics, general surgery, and more, the table features a back-lit, one-touch hand control for elevation, return-to-level, beach chair positioning, leg section movement, lateral tilt, Trendelenburg, and reverse Trendelenburg adjustments.

Specifications:

- Height: 23 in x 41 in
- Size (L × W): 79 in × 22 in (with headrest & rails)
- Table weight: 900 lb
- Weight capacity: 1200 lb lift / 1000 lb articulation
- ±30° Trendelenburg & lateral tilt
- 210° top rotation
- 40°/90° flex–reflex
- 6 in kidney bridge lift
- Back, leg, Beach Chair, and return-to-level positions
- Power Requirements: 120V, 300W

***Brochure and images are a general overview of the product. Actual product options may vary.*

VENTUREMEDICAL.COM