



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Steris Amsco 2080M Surgery Table - Refurbished



Stock Photo – Options May Vary

The Steris/Amsco 2080M manual operating table is known throughout the industry as a dependable “work horse” surgical table. The 2080M has all standard functions, $\pm$ back, $\pm$ leg, $\pm$ height, $\pm$ tilt, $\pm$ flex, $\pm$ Trendelenburg and a 300 lbs weight capacity.

The 2080M also features a manual kidney lift and manual braking system. Many accessories are available for the 2080M, including, arm boards, stirrups, anesthesia screens, shoulder positioners, urology extensions, orthopedic leg holders and extensions, and much more.

The Steris/Amsco 2080M surgical table comes completely refurbished with a comprehensive warranty and a new three piece pad set.

- Manual Operation
- Full Length Accessory Side Rails
- Manual Brakes
- $\pm$ Back, $\pm$ Leg, $\pm$ Height, $\pm$ Trendelenburg, $\pm$ Tilt, $\pm$ Flex
- Manual Kidney Lift
- 300 lbs Weight Capacity
- Includes New 3 Piece Pad
- Completely Refurbished

Specifications:

- Overall Length: 73”
- Overall Width: 21 5/8”
- Patient Surface Width: 20”
- Height Range: 27” – 45”
- Weight Capacity: 300 lbs
- Head Section Range: $\pm$ 90°
- Trend/Rev Trend: $\pm$ 20°
- Back Rest: + 55° / - 20°
- Leg Rest: + 80° / -105°
- Side Tilt: $\pm$ 10°
- Flex: 15° Back ↓ / 20° Seat ↓
- Reflex: 10° Back ↑ / 15° Seat ↑
- Pad Thickness - 2”

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReQuip, Inc.