

Steris Amsco 2080RC Surgery Table - Refurbished



Stock Photo – Options May Vary

- Wired Hand Control
- Full Length Side Rails
- Manual Brakes
- $\pm$ Back, $\pm$ Leg, $\pm$ Height, $\pm$ Trendelenburg, $\pm$ Tilt, $\pm$ Flex
- Manual Kidney Lift
- 300 lbs Weight Capacity
- Includes New 3 Piece Pad
- Completely Refurbished

The Steris/Amsco 2080RC operating table is known throughout the industry as a dependable “work horse” surgical table. The 2080RC has all standard functions, $\pm$ back, $\pm$ leg, $\pm$ height, $\pm$ tilt, $\pm$ flex, $\pm$ Trendelenburg and a 300 lbs weight capacity.

The 2080RC also features a wired hand control and manual braking system. Many accessories are available for the 2080RC, including, arm boards, stirrups, anesthesia screens, shoulder positioners, urology extensions, orthopedic leg holders and extensions, and much more.

The Steris/Amsco 2080RC surgical table comes completely refurbished with a comprehensive warranty and a new three piece pad set.

Specifications:

- Overall Length: 73”
- Overall Width: 21 5/8”
- Patient Surface Width: 20”
- Height Range: 27” – 45”
- Weight Capacity: 300 lbs
- Head Section Range: $\pm 90^\circ$
- Trend/Rev Trend: $\pm 20^\circ$
- Back Rest: $+ 55^\circ / - 20^\circ$
- Leg Rest: $+ 80^\circ / - 105^\circ$
- Side Tilt: $\pm 10^\circ$
- Flex: 15° Back $\downarrow$ / 20° Seat $\downarrow$
- Reflex: 10° Back $\uparrow$ / 15° Seat $\uparrow$
- Power: 120 Volts , 10 Amps
- Pad Thickness: 2”

***Brochure and images are a general overview of the product.
Actual product options may vary.*