

Steris/Amsco 4085 Surgery Table - Refurbished



Stock Photo – Options May Vary

The refurbished Steris 4085 Surgical Table offers a strong, reliable platform for a wide range of surgical procedures. With a 1,100 lb weight capacity in the normal position, it's built to support patients safely and comfortably in demanding clinical environments.

Key features include a sliding tabletop for C-arm imaging, wired hand control, hydro-electric brakes, and battery operation for flexibility and ease of use. The table supports multiple movements such as Trendelenburg, lateral tilt, height adjustment, powered kidney lift, and adjustable head and leg sections.

Each unit includes a new pad set and is fully reconditioned to ensure dependable performance. With its modern design, accessory compatibility, and versatile functionality, the Steris 4085 is an excellent choice for hospitals and surgical centers seeking value without compromising quality.

- Wired Hand Control
- Accessory Side Rails
- C-Arm Compatible Top
- Electro-Hydraulic Brakes
- $\pm$ Back, $\pm$ Leg, $\pm$ Height, $\pm$ flex, $\pm$ Trendelenburg, $\pm$ Tilt
- Integrated Kidney Lift
- Return-to-Level Function
- Battery Operated
- 1,100 lbs Weight Capacity (Normal Position)
- Includes New 3 Piece Pad
- Completely Refurbished

Specifications:

- Overall Length: 81"
- Overall Width: 21 $\frac{5}{8}$ "
- Patient Surface Width: 20"
- Height Range: 26" - 45"
- Weight Capacity: 1100 lbs
- Trend/Rev Trend: $\pm$ 30°
- Head Section: + 90° / - 90°
- Back Section: + 80° / - 40°
- Leg Section: 0° / -105°
- Lateral Tilt (left/right): $\pm$ 20°
- Flex/Reflex: 140°/100°
- Power: 120 Volts, 4.5 Amps
- Pad Thickness: 2"

***Brochure and images are a general overview of the product.
Actual product options may vary.*