



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Steris 5085 Surgical Table - Refurbished



Stock Photo – Options May Vary

- 18" sliding tabletop
- 1,200 lb. weight capacity (centered column)
- Large imaging area
- Low height for optimal positioning
- Removable leg section (no protruding pins)
- Telescoping stainless-steel column
- Self-leveling floor locks
- Lightweight, impact-resistant design
- LCD hand control
- Includes New Pad Set (2)

The refurbished Steris 5085 Surgical Table delivers exceptional versatility and performance for a wide range of surgical specialties. With an 18" sliding tabletop and up to 1,200 lb patient support capacity on a centered column, it eliminates the need for multiple specialty tables while ensuring safe, stable positioning.

Key features include a large imaging area for C-arm access, telescoping stainless-steel column for easier cleaning, self-leveling floor locks for stability, and an easy-to-read LCD hand control for precise adjustments. The removable leg section, low height range, and advanced articulation capabilities enhance surgical access and flexibility.

Each unit includes a new pad set and is fully refurbished to provide reliable, long-lasting performance. Designed for durability, efficiency, and clinical adaptability, the Steris 5085 is an excellent solution for hospitals and surgery centers seeking premium functionality at exceptional value.

Specifications:

- Overall Length: 81"
- Overall Width: 20"
- Height Range: 22.8" – 42.5"
- Tabletop Slide: 18" (9" head / 9" foot)
- Weight Capacity (Centered Column): 1,200 lbs
- Weight Capacity (Full Articulation): 1,000 lbs
- Weight Capacity (Articulation + Slide): 600 lbs
- Trend/Rev Trend: $\pm 30^\circ$
- Lateral Tilt (left/right): $\pm 20^\circ$
- Head Section: $+90^\circ / -90^\circ$
- Leg Section: $0^\circ / -105^\circ$

***Brochure and images are a general overview of the product. Actual product options may vary.*

VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReQuip, Inc.