

Steris CMAX Surgical Table - Refurbished



Stock Photo – Options May Vary

- C-Arm Compatible Tabletop
- Wired Hand Control
- Battery Operation
- Hydraulic Floor Locks
- Return to Level
- North American Side Rails Fit Most Accessories
- 1,100 lbs Weight Capacity in Normal Position
- Completely Refurbished
- Comprehensive Warranty

The refurbished Steris CMAX surgical table delivers strength and modern functionality at an accessible price point. With a 1,100 lb weight capacity and hydro-electric brakes, it supports a wide range of procedures across clinical environments.

Key features include wired hand control, battery operation, a sliding tabletop for C-arm compatibility, and a full range of powered movements. Movements include backrest, leg section, Trendelenburg, lateral tilt, kidney lift, and adjustable headrest. Return to Level, North American side rails, and compatibility with most surgical accessories round out the feature set.

This unit ships refurbished with a new pad set included, making it a cost-effective, multi-purpose solution ready for almost any surgical setting.

Specifications:

- Overall Length: 81"
- Overall Width: 21 5/8"
- Height Range: 26" – 45"
- Tabletop Slide: 18" (9" head / 9" foot)
- Weight Capacity (Normal Position): 1,100 lbs
- Weight Capacity (Reverse Orientation): 500 lbs
- Trend/Rev Trend: $\pm 30^\circ$
- Lateral Tilt (left/right): $\pm 20^\circ$
- Flex/Reflex: 140° / 100°
- Leg Section: 0° / -105°

***Brochure and images are a general overview of the product.
Actual product options may vary.*