



www.venturemedical.com

6008 Bonacker Drive

Tampa, FL 33610

800-627-3215

813-740-8640

Fax: 800-466-1251

sales@venturemedical.com

Steris / Hausted Surgi-Chair ESCEYE-ST - Refurbished



Stock Photo – Options May Vary

Steris Hausted Surgi-Chair ESCEYE-ST Refurbished. The Hausted Surgi Chair is a multi-position chair and stretcher intended for a patient's treatment, transport, or recovery. With the uniquely designed articulating head section, this unit can be used for specialty procedures that require precision head positioning.

The 500 lbs weight capacity can handle large patients with ease.

The Steris / Hausted SurgiChair ESC-EYE comes completely refurbished with a comprehensive warranty and a new pad.

Features:

- 5" Tente brand casters
- Four-Wheel Brake and Steer Caster System
- Thermoformed base cover
- 24 VDC Backup Battery Operation
- 500 lbs Weight Capacity
- Emergency Backrest release lever
- Includes: New Pad Set, Hand Control, New Batteries, and Manual
- Completely Refurbished

Specifications:

- Overall Length: 75"
- Overall Width: 29.4"
- Patient Surface: 22" W x 75" L
- Height Range: 22.25" – 30.5"
- Weight Capacity: 500 lbs

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM

© 2026 Venture Medical ReQuip, Inc.