



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Stryker 1001 Advantage Stretcher - Refurbished



Stock Photo – Options May Vary

The Stryker 1001 Stretcher Stretchers are perfect for emergency departments, transport, PACU, procedures, or anywhere else a versatile stretcher is needed. The 1001 has a flexible design with a 26" patient surface that can be configured to best suit your department's needs.

Standard features include crank backrest (Fowler), Trendelenburg, reverse Trendelenburg, and height adjustment. The 500 lbs weight capacity can handle large patients with ease.

The Stryker 1001 stretcher comes completely refurbished with a comprehensive warranty and a new pad.

- Gas Spring Assist Backrest
- Glideaway Side Rails
- Central Brake and Steer
- $\pm$ Height, $\pm$ Trendelenburg
- Built-in O2 Tank Holder
- 500 lbs Weight Capacity
- Includes New 4" Pad
- Completely Refurbished

Specifications:

- Overall Length: 83"
- Overall Width: 31.5"
- Patient Surface Width: 26"
- Height Range: 22" – 35.5"
- Weight Capacity: 500 lbs
- Trend/Rev Trend: $\pm 18^\circ$
- Side Rails: 14" x 54"
- Wheel Diameter: 8"
- Pad Thickness: 4"

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM