



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Stryker 1020 Trauma Stretcher - Refurbished



Stock Photo – Options May Vary

The Stryker 1020 Trauma Stretcher offers stability for the critically injured. Stryker's trauma stretcher provides a secure platform for assessing and tending to critically injured and immobile patients. Its innovative design allows for full-length X-ray capability without having to transfer or reposition the patient.

The 1020 Stretcher features full-length x-ray capability, securely-mounted, lateral c-spine cassette holder, and precisely position slide-out trays for supreme imaging results as well as a 90-degree backrest that helps facilitate imaging access.

The Stryker 1020 Trauma stretcher comes completely refurbished with a comprehensive warranty and a new pad.

- Full-length X-ray Cassette Tray
- Dual Cylinder Pneumatic Backrest
- 500 lbs Weight Capacity
- Fold-down Siderails
- $\pm$ Height, $\pm$ Trendelenburg
- Four IV receptacles
- Full-length X-ray Cassette Tray
- Integrated Storage Tray
- Built-in O2 Tank Holder
- Roller Bumpers
- X-ray Centering Label
- Retractable Fifth-wheel Steering
- Includes New 4" Pad
- Completely Refurbished

Specifications:

- Overall Length: 83"
- Overall Width: 31.5"
- Patient Surface Width: 26" W x 76"L
- Weight Capacity: 500 lbs
- Height Range: 38.5" - 24.75"
- Litter Positioning: Backrest 0° - 90°
- Trend/Reverse Trend. $\pm$ 18°
- Siderails: 13"H x 55"L
- Caster Diameter: 8"
- Pad Thickness: 4"

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReQuip, Inc.