



www.venturemedical.com

Venture Medical ReQuip, Inc.  
6008 Bonacker Drive  
Tampa, FL 33610  
800-627-3215 | 813-740-8640  
Fax: 800-466-1251  
sales@venturemedical.com

## Stryker Gynnie 1061 OB-GYN Stretcher - Refurbished



*Stock Photo – Options May Vary*

- Pneumatic Assist Fowler
- Glide-Away Side Rails
- Central Brake and Steer
- $\pm$  Height,  $\pm$  Trendelenburg
- 500 lbs Weight Capacity
- Built-in O2 Tank Holder
- Includes New 4" Pad
- Built-in Stirrups
- Completely Refurbished

The Stryker Gynnie 1061 OB-GYN Stretcher is perfect for emergency departments, labor and delivery, or gynecological procedures. The Stryker 1061 features a pneumatic assist fowler (backrest), folding leg section, and built-in stirrups.

Additional functions include Trendelenburg, reverse Trendelenburg, along with central brake and steer, glide-away side rails, and integrated O2 tank well.

The Stryker Gynnie 1061 OB-GYN stretcher comes completely refurbished with a comprehensive warranty and a new pad.

### Specifications:

- Overall Length: 81½"
- Overall Width: 31½"
- Patient Surface Width: 26"
- Height Range: 22.5" – 36"
- Weight Capacity: 500 lbs
- Fowler Angle: 0° - 90°
- Trend/Rev Trend:  $\pm$  18°
- Side Rails: 13" x 55"
- Wheel Diameter: 8"
- Pad Thickness: 4"

*\*\*Brochure and images are a general overview of the product.  
Actual product options may vary.*

**VENTUREMEDICAL.COM**

© 2007-2026 Venture Medical ReQuip, Inc.