



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Stryker 1080 C-Arm Stretcher - Refurbished



Stock Photo – Options May Vary

- Full-Length C-Arm Access
- Fold-Down Side Rails
- Central Brake and Steer
- $\pm$ Height, $\pm$ Trendelenburg
- Pneumatic Assist Fowler (x 2)
- End Foot Controls
- 500 lbs Weight Capacity
- Includes New 2" Pad
- Completely Refurbished

The Stryker 1080 Fluoroscopy Stretcher is perfect for patients and procedures that require C-Arm imaging.

The Stryker 1080 features dual pneumatic assist backrests, a U-shaped base for ease of access, and a radiolucent litter top.

Additional features include Trendelenburg and Reverse Trendelenburg along with Central Brake/Steering, Fold-Down Side Rails, and end foot controls.

The Stryker 1080 Fluoroscopy stretcher comes completely refurbished with a comprehensive warranty and a new pad.

Specifications:

- Overall Length: 87"
- Overall Width: 30.5"
- Weight Capacity: 500 lbs
- Patient Surface Width: 25.5"
- Fluoro Window: 19.5" – 65.5"
- Height Range: 24" – 33.5"
- Fowler Angle: 0° - 90°
- Trend/Rev Trend: $\pm$ 8°
- Side Rails: 13" x 55"
- Wheel Diameter: 6"
- Pad Thickness: 2"

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReQuip, Inc.