



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Stryker 1105 Prime Stretcher - Refurbished



Stock Photo – Options May Vary

The Stryker 1105 Prime Series Stretchers are extraordinarily mobile and easy to use. 5th-Wheel steering helps make transport safe and efficient while reducing the physical strain on the clinician.

Features include Dual-side mounted foot controls with uni-lower pedal, dual end brake/steer controls, glideaway siderails with integrated push/pull handles, and pneumatic backrest.

The Stryker 1105 Prime stretcher comes completely refurbished with a comprehensive warranty and a new pad.

- Glideaway Side Rails
- Central Brake and Steer
- $\pm$ Height, $\pm$ Trendelenburg
- Pneumatic Backrest
- O2 Tank Storage
- 700 lbs Weight Capacity
- IV Pole (optional)
- Includes New 4" Pad
- Completely Refurbished

Specifications:

- Overall Length: 85"
- Overall Width: 34"
- Patient Surface Width: 26"
- Height Range: 20.75"– 34"
- Weight Capacity: 700 lbs
- Fowler Angle: 0° - 90°
- Trend/Rev: +16° to -17°
- Wheel Diameter: 8"
- Pad Thickness: 4"

***Brochure and images are a general overview of the product. Actual product options may vary.*

VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReQuip, Inc.