



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Stryker 5050 Stretcher Chair - Refurbished



Stock Photo – Options May Vary

The Stryker 5050 combines the advantages of a stretcher with that of a recovery recliner. This versatile stretcher chair features a 400 lbs weight capacity, central brake and steer function, O2 tank well and fold down side rails with arm pads. Dual pedestal hydraulics increase stability and the footrest adjusts to various heights to accommodate each individual patient.

Functions include: hydraulic height adjustment, trendelenburg and reverse trendelenburg and a gas spring back rest adjustment.

The Stryker 5050 stretcher chair comes completely refurbished with a comprehensive warranty and a new pad.

- Gas Spring Assist Backrest
- Fold Down Side Rails
- Central Brake and Steer
- $\pm$ Height, $\pm$ Trendelenburg
- Built-in O2 Tank Holder
- 400 lbs Weight Capacity
- Includes New 3" Pad
- Completely Refurbished

Specifications:

- Overall Length: 76"
- Overall Width: 30"
- Patient Surface Width: 24"
- Height Range: 22" - 33.5"
- Weight Capacity: 400 lbs
- Trend/Rev Trend: $\pm 18^\circ$
- Backrest Adjustment: 0° - 90°
- Side Rails: 10" x 31"
- Wheel Diameter: 6"
- Pad Thickness: 3"

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReQuip, Inc.