



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Stryker 5051 Eye Stretcher Chair - Refurbished



Stock Photo – Options May Vary

The Stryker 5051 Eye Stretcher Chair is ideal for surgery centers and small spaces. Stryker's eye stretcher chairs maximize patient comfort throughout procedures without compromising patient care and outcomes.

The 5051 also features all the benefits of the standard stretcher chair, plus dual-articulation, enhanced-clearance headpiece, adjustable height footrest, 24" wide patient surface, and a fixed position push bar.

Standard features include: dual cylinder pneumatic backrest and dual pedestal hydraulics, swing-down siderails with comfort armrest, and thermoformed ABD hood with integrated storage tray.

The Stryker 5051 Eye Stretcher comes completely refurbished with a comprehensive warranty and a new pad.

- Dual-Articulation Headpiece
- Head Pads - Ultra Comfort
- Includes New 2" Vinyl Mattress
- 6" Casters
- 24" Wide Patient Surface
- Adjustable Height Footrest
- Dual Cylinder Pneumatic Backrest
- Dual Pedestal Hydraulics
- Fixed Position Push Bar
- Four IV Receptacles
- Hands-free Trendelenburg
- Completely Refurbished

Specifications:

- Overall Dimensions: 76"L x 30"W
- Litter Width: 27.5"
- Weight Capacity: 400 lbs
- Height Range: 22" to 33.5"
- Litter Positioning:
 - Backrest: 0° to 90°
 - Leg: 0° to 80°
- Trend./Reverse Trendelenburg : $\pm 18^\circ$
- Head Pads:
 - Flat: 3", 2", or 1"
 - Concave: 3" or 2"
- Siderails:
 - Standard: 10" to 31"
 - Transfer: 9" to 31"
- Patient Surface: 24" x 74"
- Caster Diameter: 6"

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReQuip, Inc.