

Stryker Atlas 660 Stretcher - Refurbished



Stock Photo – Options May Vary

The Stryker Atlas 660 Stretcher is the gurney solution for all departments. The 660-pound weight capacity with an extra-wide 31" patient top meets the needs of all patients. The Stryker Big Wheel technology is proven to offer easier mobility and control. It requires 50% less maximum startup force and 60% less turning effort.

Additional features include a power-washable design with 660 lbs weight capacity, 4 integrated IV receptacles, a litter-mounted O2 holder, and glide-away siderails.

The Atlas 660 features all standard features plus big wheel steering to provide superior maneuverability and cornering.

The Stryker Atlas 660 Stretcher comes completely refurbished with a comprehensive warranty and a new pad.

- Gas Spring Assist Backrest
- Glide-away Side Rails
- Central Brake and Steer
- $\pm$ Height
- $\pm$ Trendelenburg
- Built-in O2 Tank Holder
- 660 lbs Weight Capacity
- Includes New 4" Pad
- Completely Refurbished

Specifications:

- Overall Length: 84"
- Overall Width Siderails Up: 34.5"
- Overall Width Siderails Down: 31"
- Patient Surface Width: 31 1/2" x 76"
- Weight capacity: 650 lbs
- Height range: 22" - 33"
- Knee Gatch: 30°
- Trend/Rev Trend: $\pm$ 15°
- Caster Diameter: 8"
- Pad Thickness: 4"

***Brochure and images are a general overview of the product.
Actual product options may vary.*