



www.venturemedical.com

Venture Medical ReQuip, Inc.  
6008 Bonacker Drive  
Tampa, FL 33610  
800-627-3215 | 813-740-8640  
Fax: 800-466-1251  
sales@venturemedical.com

## Stryker Renaissance/Advantage Stretcher - Refurbished



Stryker 1731 Pictured\*

Stock Photo – Options May Vary

- Gas Spring Assist Backrest
- Glide-Away or Collapsible Side Rails
- Central Brake and Steer
- $\pm$  Height,  $\pm$  Trendelenburg
- Built-in O2 Tank Holder
- 500 lbs Weight Capacity
- Includes New 4" Pad
- Completely Refurbished

The Stryker Renaissance/Advantage Series Stretchers are perfect for emergency departments, transport, PACU, procedures, or anywhere else a versatile stretcher is needed.

Standard features include lift assist backrest (Fowler), Trendelenburg, reverse Trendelenburg, and height adjustment. The 500 lbs weight capacity can handle large patients with ease.

The Stryker Renaissance/Advantage series stretcher comes completely refurbished with a comprehensive warranty and a new pad.

### Specifications\*:

- Overall Length: 83"
- Overall Width: 31.5"
- Patient Surface Width: 26" or 29"
- Height Range: 22" - 35½"
- Weight Capacity: 500 lbs
- Trend/Rev Trend:  $\pm$  18°
- Side Rails: 14" x 54"
- Caster Size: 8"
- Pad Thickness: 4"

\*Specifications are approximate.

\*Specifications will vary by model and features selected.

\*\*Brochure and images are a general overview of the product.  
Actual product options may vary.

# VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReQuip, Inc.