



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Stryker SM204 (1015) Stretcher - Refurbished



Stock Photo – Options May Vary

The Stryker SM204 Big Wheel stretcher features mobility, simplicity, and versatility. Fully equipped and completely adaptable, the M-Series allows you to choose the mobility solution that best fits your needs. Supremely mobile and user-friendly, these stretchers greatly reduce the physical strain on clinicians without sacrificing patient care and comfort.

Additional features include: Big Wheel steering, power-washable design with 700 lbs weight capacity, 4 integrated IV receptacles, litter-mounted O2 holder, and glide-away siderails.

The SM204 features all standard features plus Big Wheel steering, which provides superior maneuverability and cornering.

The Stryker SM204 stretcher comes completely refurbished with a comprehensive warranty and a new pad.

- Gas Spring Assist Backrest
- Glide-away Side Rails
- Central Brake and Steer
- $\pm$ Height
- $\pm$ Trendelenburg
- Built-in O2 Tank Holder
- 700 lbs Weight Capacity
- Includes New 4" Pad
- Completely Refurbished

Specifications:

- Overall Length: 84"
- Overall Width: 33.5"
- Patient Surface Width: 26" x 75.25"
- Weight capacity: 700 lbs
- Height range: 20.75" - 34.5"
- Knee Gatch: 40°
- Trend/Rev Trend: $\pm$ 16°
- Side Rails: 14" x 57.75"
- Caster Diameter: 8"
- Pad Thickness: 4"

***Brochure and images are a general overview of the product. Actual product options may vary.*

VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReQuip, Inc.