



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Stryker Trio Mobile Surgical Platform/Stretcher - Refurbished



Stock Photo – Options May Vary

- All-In-One Efficiency
- Manual Operation
- Glide-away Side Rails
- Central Brake and Steer
- $\pm$ Height
- $\pm$ Backrest
- $\pm$ Trendelenburg
- $\pm$ Lateral Tilt
- 500 lbs Weight Capacity
- Includes New 3" Pad
- Completely Refurbished

The Stryker Trio (1033) mobile surgery platform virtually eliminates the need for patient transfers by allowing the patient to stay on the surface for all phases – Pre-Op, Surgery, and Post-Op. The Stryker Trio allows the user to move the unit in the common surgical features, which make it ideal for the ASC and in-office procedure environment. The Trio has manual operations, allowing for simplicity and longevity, as well as a secure braking system, including a steer function for when transporting the patient. mobility, simplicity, and versatility.

Additional features include standard accessory siderails, glide-away patient siderails, IV pole receptacles, 4-wheel central brake system, 5th wheel steering, and low-profile base. Optional accessories include drain pan, siderail pads, foot extender, headpiece extender, Eye surgery backrest, and more.

The Stryker 1033 Trio surgical table/stretcher comes completely refurbished with a comprehensive warranty and a new pad.

Specifications:

- Overall Length: 80"
- Overall Width: 33"
- Patient Surface Width: 27" x 78"
- Weight capacity: 500 lbs
- Height range: 26" - 40"
- Backrest: 1° to 80°
- Lateral Tilt: $\pm$ 15°
- Foot Section: 0° - 50°
- Trend/Rev Trend: +20°/-15°
- Side Rails: 12" x 49"
- Caster Diameter: 6"

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReQuip, Inc.