



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Wy'East Totalift II Stretcher Chair - Refurbished



Stock Photo – Options May Vary

**** Please Note: Upper Side Rails and Contoured Headrest are Optional**

- Transfer Slide Top
- Gas Spring Assist Backrest
- $\pm$ Height, + Trendelenburg
- Tuck-Away Side Rails
- 4 Wheel Locking
- Built-in O2 Tank Holder
- 400 lbs Weight Capacity
- Includes New 2" Pad
- Completely Refurbished

The Wy'East Totalift II Stretcher Chair is an extremely versatile stretcher that is perfect for a wide variety of departments within your medical facility.

The Totalift Stretcher features a pneumatic assist fowler (backrest), hydraulic height adjustment, central brake and steering system, side rails, O2 tank well, and sliding patient surface for easy transfers.

The Wy'East Totalift II comes completely refurbished with a comprehensive warranty and a new pad.

Specifications:

- Overall Length: 71½"
- Overall Width: 27"
- Height Range: 23" - 36"
- Weight Capacity: 400 lbs
- Fowler Angle: 0° - 80°
- Trendelenburg: 20°
- Wheel Diameter: 4"
- Pad Thickness: 2"

***Brochure and images are a general overview of the product. Actual product options may vary.*

VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReQuip, Inc.