

ATS[®]
2000
ATS

THE ATS 2000
TOURNIQUET SYSTEM



*The Choice for
Safety, Reliability
and Convenience*



Safety.

Zimmer's dual-port system
helps protect your patients by
monitoring and controlling actual
cuff/ bladder pressures.

Reliability.

Advanced microprocessor technology
includes self-diagnostic tests and
continuous monitoring for consistently
dependable performance.

Convenience.

Fast inflation and deflation save time
for busy operating room personnel.



The safety The reliability and co

Dual controls, dual cuffs, and dual readouts provide independent and precise monitoring and performance for bilateral or Bier Block procedures.

Accurate pressure readings are taken directly from each cuff.



Color-coded cuffs from 8 in. to 42 in. accommodate virtually all patients and limbs.



Easy-to-read display adjusts easily, too.

Simple setup includes user-defined default pressure and time settings.



you need.

convenience you want.



Low-profile cuffs stay in place without discomfort.

Select reusable cuffs or disposable sterile cuffs to eliminate cross-contamination.

Cuff sleeves help reduce wrinkles, pinching, and shearing of soft tissue.



Positive-locking connectors supply safer, leak-free inflation.

Angled ports keep connectors clear of the surgical site.

Visual and adjustable audio alarms readily detect pressure variations.



Lightweight portable unit can be mounted on an I.V. pole or placed on a table.

One-hour battery backup – just in case.



The A.T.S. 2000 Tourniquet System

Combining safety, reliability and convenience so you can focus on surgery.

The built-in safety, reliability and convenience of the A.T.S. 2000 Tourniquet frees you to focus on surgery. Dual displays that monitor actual cuff pressure, precise microprocessor controls, and one-hour battery backup help assure patient safety. User-set alarms warn of unacceptable or hazardous conditions.

The dual-cuff display is easy to read, even from a distance, and touch pad controls are simple to read and use. Fast cuff inflation and deflation save time. A large selection of cuffs, disposable and reusable, fits virtually all patients. The lightweight unit can be conveniently positioned, and positive-locking connectors have angled ports to keep lines clear of the surgical site.



Ordering Information

A.T.S. 2000 TOURNIQUET SYSTEM

Prod. No.	Description
60-2000-101	A.T.S. 2000 With Hoses Only. No Stand, Basket or Cuffs.
60-2000-102	A.T.S. 2000 With Hoses and A.T.S. Cylindrical Cuffs (60-7500-001/007).
60-2000-103	A.T.S. 2000 With Hoses, A.T.S. Cylindrical Cuffs (60-7500-001/007), Stand and Basket.
60-2000-104	A.T.S. 2000 With Hoses, Stand and Basket. No Cuffs.

STERILE DUAL HOSE EXTENSION WITH PLC (BOX 20)

Prod. No.	Size
60-1812-101	36 in. (91cm)

DISPOSABLE STERILE A.T.S. CYLINDRICAL CUFFS

DUAL PORT, SINGLE BLADDER (BOX 10)

Without Sleeve	With Sleeve	Size
60-7070-101	60-7070-151	8 in. (20cm) x 2.75 in. (7cm)
60-7070-102	60-7070-152	12 in. (30cm) x 3.5 in. (9cm)
60-7070-103	60-7070-153	18 in. (46cm) x 4 in. (10cm)
60-7070-104	60-7070-154	24 in. (61cm) x 4 in. (10cm)
60-7070-105	60-7070-155	30 in. (76cm) x 4 in. (10cm)
60-7070-106	60-7070-156	34 in. (86cm) x 4 in. (10cm)
60-7070-107	60-7070-157	42 in. (107cm) x 4 in. (10cm)

DUAL PORT, SINGLE BLADDER (BOX 1)

Without Sleeve	With Sleeve	Size
-----	60-7070-153-01	18 in. (46cm) x 4 in. (10cm)
-----	60-7070-154-01	24 in. (61cm) x 4 in. (10cm)
-----	60-7070-155-01	30 in. (76cm) x 4 in. (10cm)
-----	60-7070-156-01	34 in. (86cm) x 4 in. (10cm)

DUAL PORT, DUAL BLADDER (BOX 10)

Without Sleeve	With Sleeve	Size
60-7080-101	60-7080-151	12 in. (30cm) x 3.5 in. (9cm)
60-7080-102	60-7080-152	18 in. (46cm) x 4 in. (10cm)
60-7080-103	60-7080-153	24 in. (61cm) x 4 in. (10cm)

REUSABLE A.T.S. CYLINDRICAL CUFFS

DUAL PORT, SINGLE BLADDER (BOX 1)

Prod. No.	Size
60-7500-001	8 in. (20cm) x 2.75 in. (7cm)
60-7500-002	12 in. (30cm) x 3.5 in. (9cm)
60-7500-003	18 in. (46cm) x 4 in. (10cm)
60-7500-004	24 in. (61cm) x 4 in. (10cm)
60-7500-005	30 in. (76cm) x 4 in. (10cm)
60-7500-006	34 in. (86cm) x 4 in. (10cm)
60-7500-007	42 in. (107cm) x 4 in. (10cm)

DUAL PORT, DUAL BLADDER (BOX 1)

Prod. No.	Size
60-7555-001	12 in. (30cm) x 3.5 in. (9cm)
60-7555-002	18 in. (46cm) x 4 in. (10cm)
60-7555-003	24 in. (61cm) x 4 in. (10cm)

REUSABLE CONTOUR CUFFS

DUAL PORT, SINGLE BLADDER (BOX 1)

Prod. No.	Description	Width
60-7450-002	Upper Limb Cuff	5.75 in. (15cm)
60-7450-003	Thigh Cuff	5.75 in. (15cm)
60-7450-004	Lower Leg Cuff	5.75 in. (15cm)

DUAL PORT, DUAL BLADDER (BOX 1)

Prod. No.	Description	Width
60-7450-011	Upper Limb IVRA Cuff	5.75 in. (15cm)

REUSABLE SMALL ADULT CUFFS

DUAL PORT, SINGLE BLADDER (BOX 1)

Prod. No.	Description	Size
60-8100-001	XSm Adult Cuff	15 in. (38cm) x 3 in. (8cm)
60-8100-002	Sm Adult Cuff	17.5 in. (44cm) x 3.5 in. (9cm)

CUFF SLEEVES ONLY (BOX 10)*

FOR A.T.S. REUSABLE AND DISPOSABLE CYLINDRICAL CUFFS

Prod. No.	Description
60-8000-008	8 in. (20cm) Cuff Sleeve
60-8000-012	12 in. (30cm) Cuff Sleeve
60-8000-018	18 in. (46cm) Cuff Sleeve
60-8000-024	24 in. (61cm) Cuff Sleeve
60-8000-030	30 in. (76cm) Cuff Sleeve
60-8000-034	34 in. (86cm) Cuff Sleeve
60-8000-042	42 in. (107cm) Cuff Sleeve

*Sold Sterile

ACCESSORIES

Prod. No.	Description
60-4017-001	Dual Cuff Hose Blue/White
60-4018-001	Dual Cuff Hose Red/White
60-4022-001	Tourniquet Stand (Includes Base and Pole)
60-1908-001	Accessory Basket

Above cuffs come standard with Positive Locking Connectors (PLC). Cuffs with Luer Lock Connectors also available.