

Worth closer examination



Dual Luminance
Settings



Mobile Base MI-300



Brilliant, Color-
Accurate Illumination

MI-300

LED Examination Lights

Presenting the best in high-value, space-saving, color-accurate LED exam lights. The MI-300 offers two brightness settings and its mobile base allows it to be moved anywhere easily. Both exam lights provide superior brightness and clarity for more exacting diagnoses.

Specifications

Light Intensity	
MI-300.....	30,000 or 15,000 lux at 24" (60 cm) 50,000 lux or 25,000 lux at 18" (45.7 cm)
Color Temperature.....	4,300° K
CRI.....	.97
LED Type/Life.....	CoB / 50,000 hours
Operating Voltage.....	100VAC-240VAC/50-60Hz
Light Head Input Voltage.....	24VDC
Unit Height.....	MI-300: 66.0" (1,625mm) MI-150: 64.0" (1,676mm)
Warranty.....	3 years (limited)



19749 Dearborn Street, Chatsworth, CA 91311
(800) 831-1222 | www.medillum.com