



**COVIDIEN**

*positive results for life™*

## RapidVac™ Smoke Evacuator System

The RapidVac™ smoke evacuator system improves your environment, visibility\* and performance — *quietly*.

### **MORE POWER**

The RapidVac™ smoke evacuator system has 30 percent more power at maximum air flow than our previous model.

### **EFFECTIVE FILTER TECHNOLOGY**

All-in-one, four-stage, high-technology, long-life, three-port, RFID-equipped ULPA filter provides 99.9995 percent filtration efficiency at 0.12 micron particle size.

### **EASY TO USE\***

Easily understandable graphics, easy to set up, easily determined airflow and filter life, and easy filter replacement.

### **ECONOMICAL**

Long filter life, 25 hours activation time. Three ports for common filter sizes. No need for prefilters or adapters, reducing cost per procedure.

### **QUIET\***

The RapidVac™ smoke evacuator system is quieter than our previous model and judged in surveys as quieter than expected in the O.R.



### **RAPIDVAC™ SMOKE EVACUATOR SYSTEM**

The RapidVac™ smoke evacuator teams up with the ForceTriad™ energy platform and the Force FX™ and Force EZ™ generators. Optional accessories are available to enable other generators to synchronize with RapidVac™.

## RapidVac™ Smoke Evacuator Technical Specifications

<b>Dimensions:</b>	8 in. high x 14 in. wide x 18 ¾ in. deep
<b>Weight:</b>	20 lbs.
<b>Filter:</b>	ULPA: 99.9995% efficient at 0.12 microns Pre-filter: Rated for 0.3 micron filtration
<b>Electrical:</b>	Voltage operating range: 200-132 VAC, 50 Hz; current: 10 amps, maximum 200-240 VAC, 60 Hz; current: 8 amps, maximum
<b>Compatible Carts:</b>	FT900 Mounting cart system - main, upper, or lower shelf UC8009 Universal mounting cart UC8013 Suspended shelf for
<b>Certifications and Markings:</b>	UL, CE, cUL

### ORDER INFORMATION

Catalog Number	Description	Order Quantity
<b>Smoke Evacuator and Accessories</b>		
SE3690	RapidVac™ smoke evacuator 100-132 V	1 each
SE3695	RapidVac™ smoke evacuator 200-240 V	1 each
SEA3730	RapidVac™ system interlink cable, for use with the ForceTriad™	1 each
SEA3740	RapidVac™ remote switch activator (for use with any electrosurgical generator)	1 each
SEA3745	RapidVac™ footswitch, pneumatic, with 3 m (10 ft) cord	1 each
<b>Tubeing and Filters</b>		
SEA3700	RapidVac™ smoke evacuator ULPA filters (2 pack)	2/case
SEA3735	AccuVac™ pencil adapter attachment, with 1 cm (3/8 in.) tubing, 3 m (10 ft), sterile	Coming soon
SEA3720	Smoke vac tubing kit w/valve, sterile	12/case
SEA3715	2.2 cm (7/8 in.) tubing (10 ft), sterile	12/case
SEA3710	1 cm (3/8 in.) tubing, 3 m (10 ft), sterile	10/case
SEA3725	Laser-resistant wand, 2.2 cm (7/8 in.), sterile	Coming soon
<b>AccuVac™ Smoke Evacuation Pencil Kits</b>		
E2350HS**	Rocker switch pencil with EDGE™ coated blade, smoke evacuation attachment, and holster	25/case
E2450HS**	Button switch pencil with EDGE™ coated blade, smoke evacuation attachment, and holster	25/case
E2515HS**	Rocker switch pencil with smoke evacuation, attachment and holster	25/case
E2516HS**	Button switch pencil with smoke evacuation, attachment and holster	25/case
E3590	AccuVac™ smoke evacuation attachment	25/case
E3595	AccuVac™ smoke evacuation attachment extender	10/case

\*In direct postsurgical surveys of surgeons and nurses on their use of the RapidVac™ Smoke Evacuator, the majority affirmed support for the following: (1) Easy to set up and easy to connect tubing (32 out of 32); (2) Filter easy to install (32 out of 32); (3) Filter easy to remove (31 out of 31); (4) Easy to adjust outflow (31 out of 31); (5) Quieter in the O.R. than expected (20 out of 20); (6) Improved visibility (31 out of 31). Data on file at Covidien.

\*Requires SEA3710

\*\*Requires SEA3710



COVIDIEN, COVIDIEN with logo and "positive results for life" are trademarks of Covidien AG. © 2009 Covidien. All rights reserved.

R0007519 Rev. 2009/04