

# GLIDESCOPE®

## AVL Single-Use

When time is critical, make your go-to device the GlideScope® AVL.

Ideal for demanding OR, ED, and ICU settings, the AVL provides a consistently clear view of the airway and enhanced maneuverability, enabling fast intubations.

The AVL Single-Use system gives you 6 sizes to cover a wide range of patients.

Add real-time recording and an onboard tutorial, and you'll see why GlideScope AVL is the one you go to for advanced video laryngoscopy.




### AVL=Advanced Airway Views

- ▶ Large digital color monitor
- ▶ Innovative blade sizing, angulation, and camera positioning
- ▶ Reveal™ anti-fog feature, with a rapid heating profile to resist lens fogging

### AVL=Advanced User Features

- ▶ Real-time snapshot and video recording directly to a USB flash drive helps confirm tube placement and facilitate teaching
- ▶ Onboard *GlideScope 4-Step Technique* tutorial makes AVL easy to use, learn, and teach
- ▶ Ergonomic design and rugged construction for multiple clinical applications
- ▶ 6 sterile, single-use Stat sizes available, designed for a broad range of patient sizes



## Applications

- First-use intubations, replacing direct laryngoscopy
- Normal or restricted oropharyngeal views
- Cormack-Lehane Grades I–IV laryngeal views
- Trauma airways—excellent when dealing with blood and secretions in the airway
- Nasal tracheal intubation
- Video-guided tube exchange in the ICU

## System Components

- Digital, large color monitor
- Reusable video baton with digital camera
- Single-use Stat, choice of 6 sizes (box of 10)
- GlideRite® Rigid Stylet
- Protective video baton cradle
- Choice of mobile stand or IV pole mount
- Operations & Maintenance Manual
- Instructional In-service video
- 1-year warranty

## GlideScope® Products and Systems

| SINGLE-USE SYSTEM                                                                   |                                                                                     |                                                                                       |                                                                                       |                                                                                      |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| VIDEO BATON 1-2<br>For GVL® 0 - 2.5 Stats                                           | GVL 0                                                                               | GVL 1                                                                                 | GVL 2                                                                                 | GVL 2.5                                                                              |
|   |   |     |   |  |
| VIDEO BATON 3-4<br>For GVL 3 - 4 Stats                                              | GVL 3                                                                               |                                                                                       | GVL 4                                                                                 |                                                                                      |
|  |  |                                                                                       |  |                                                                                      |
| GLIDESCOPE PREMIUM CART                                                             |                                                                                     | AVL PORTABLE STAND                                                                    |                                                                                       |                                                                                      |
|  |                                                                                     |  |                                                                                       |                                                                                      |

Note: System configurations vary.



Verathon Inc.  
20001 North Creek Parkway  
Bothell, WA 98011  
USA  
Tel: +1 800 331 2313  
(USA & Canada only)  
Tel: +1 425 867 1348  
Fax: +1 425 883 2896



Verathon Medical (Europe) B.V.  
Willem Fenengastrat 13  
1096 BL Amsterdam  
The Netherlands  
Tel: +31 (0) 20 210 30 91  
Fax: +31 (0) 20 210 30 92

Verathon Medical (Australia) Pty Limited  
Unit 9, 39 Herbert Street  
St Leonards NSW 2065  
Australia  
Tel: 1800 613 603 (Within Australia)  
Fax: 1800 657 970  
Tel: +61 2 9431 2000 (International)  
Fax: +61 2 9475 1201