

AcralMed Stand Alone Leg Prepper | OTA154



Stock Photo - Options May Vary

FACTORY NEW

AcralMed Stand Alone Leg Prepper | OTA154 (also known as MCM154) is a leg prep stand that supports a patient's leg during skin prep without attaching to the table. Most leg preppers mount to the table rail, which limits positioning. This one works beside the table or away from it.

It also works on stretchers and gurneys, which have no rail. Meanwhile, the base can slide under the OR table, and the upper portion swivels 360 degrees to reach the right position. A thumb knob adjusts the height.

As a result, you can place the stand beside the table, off the end, or wherever you need it.

FEATURES

- Stand alone design, with no table rail needed
- Works with any OR table, stretcher, or gurney
- Upper portion swivels 360 degrees
- Thumb knob adjusts the height
- Base slides under the OR table

SPECIFICATIONS

- SKU: OTA154 (also MCM154)
- Rotation: Upper portion swivels 360 degrees
- Height Adjustment: Thumb knob
- Base: Slides under the OR table
- Use With: Any OR table, stretcher, or gurney

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)